

BARIATRIC ECONOMY BED SYSTEMS



(Shown with optional side rails and mattress)

Model 3150 BARI LOW EC BED Weight Capacity: 850 Pounds

Rigorously tested: ZERO Failures

| FEATURE | BENEFIT |
|---|--|
| <ul style="list-style-type: none"> ● 9.75" low deck to floor height | <ul style="list-style-type: none"> Minimizes injuries due to falls and provides easy patient access and egress |
| <ul style="list-style-type: none"> ● 26" high deck position | <ul style="list-style-type: none"> Comfortable height for caregiver |
| <ul style="list-style-type: none"> ● Three function pendant control | <ul style="list-style-type: none"> Convenient clinical positioning and patient comfort including fowler positioning, leg lift and deck high/low |
| <ul style="list-style-type: none"> ● 84" deck length | <ul style="list-style-type: none"> For your taller patients |
| <ul style="list-style-type: none"> ● 36"-42"-48" adjustable deck width | <ul style="list-style-type: none"> Narrows for transport, wide enough for most bariatric patients |
| <ul style="list-style-type: none"> ● Two-piece frame construction | <ul style="list-style-type: none"> Enables easy delivery and setup |
| <ul style="list-style-type: none"> ● Optional foot end side rails | <ul style="list-style-type: none"> Added safety and security |
| <ul style="list-style-type: none"> ● Optional emergency battery backup | <ul style="list-style-type: none"> Provides power when A/C power is not available |

TOTAL SOLUTIONS FOR BARIATRIC PATIENT CARE