

GENDRON, INC.  
520 W Mulberry St  
Suite 100  
Bryan Ohio 43506

Effective January 1 2009

**CUSTOMER SERVICE**

Phone: 800-537-2521  
Fax: 419-636-9261  
www.gendroninc.com

**Maxi Rest Bariatric Beds  
Extra Care Bariatric Bed**

**Weight Capacity: Up To 1000 Pounds**

Date: .....	P.O. #: .....
Dealer #: .....	Contact: .....
Bill To: .....	Ship To: .....
Address: .....	Address: .....
City/State/Zip: .....	City/State/Zip: .....
Phone No: (.....).....	Phone No: (.....).....

To order, check boxes on the price list below. Add up prices, insert "Order total." Mail, fax or call-in order.

**MODEL NUMBER**

- 4054SD 3 function electric, two piece frame, head end rails, head & foot boards, width and length adjustable 1000 lb. cap. \$10200
- 4054SDX 5 function electric, two piece frame, head end rails, head & foot boards, width and length adjustable 1000 lb. cap. \$11000
- GEN 1000\* 5 function electric, one piece frame, head & foot rails, head & foot boards, width and length adjust. 1000 lb. cap. \$18300

**Note:** three function electric features height adjustment, fowler back rest and leg lift  
five function electric features height adjustment, fowler back rest, leg lift, trendelenburg and reverse trendelenburg  
\*five function electric features height adjustment, fowler back rest, drop foot end section, trendelenburg and reverse trendelenburg  
**Important:** all models adjust from 39" to 48" to 54" wide utilizing integral extension panels, length adjustable to 88"

**OPTIONAL ACCESSORIES**

**Mattress**

- 3980F 6" bariatric mattress, 39 x 80, 1100 lb ..... \$820
- 4880F 6" bariatric mattress, 48 x 80, 1100 lb. .... \$890
- 4880-1F 6" bariatric mattress, 39-48 x 80, 1100 lb. .. \$1150
- 5480-1F 6" bariatric mattress, 39-54 x 80, 1100 lb. .. \$1350
- 4148F mattress for GEN1000 39-48 x 80, 1000 lb.. \$1500
- 4154F mattress for GEN1000 39-54 x 80, 1000 lb.. \$1700

**Note:** models 4880-1F and 5480-1F mattress includes side bolsters to expand to 48" or 54" wide

**Bolsters / Foot End**

- 3908F 6"x8"x39" mattress bolster & deck panel 1100 lb. \$350
- 4808F 6"x8"x48" mattress bolster & deck panel 1100 lb. \$375
- 5408F 6"x8"x54" mattress bolster & deck panel 1100 lb. \$425
- 22155 deck extension panel only ..... \$250

**Side Rails**

- head end side rails, pair, included with bed ..... N/C.
- 4043 one pair side rails, foot end ..... \$565
- 4044-2 one pair side rails extra tall for Maxi Rest ..... \$215  
Note: foot end side rails are standard on model GEN1000

**Trapeze Frame**

- 4045 over bed trapeze and frame 80" long ..... \$1250
- 4046 over bed trapeze and frame 88" long ..... \$1450

**Scale System**

- W integral weighing scales, caster stem mount ..... \$5000  
Note: integral scales are standard on model GEN1000

**Caster Wheels**

- hospital casters, two locking ..... Standard

**CPR Release**

- CPR fowler back rest release ..... \$375  
**Note:** CPR fowler is standard on model GEN1000  
CPR release not available on split frame beds

**Accessory Wall Bumper**

- 4047 head or foot full width wall bumper SD/SDX beds...\$215

**STANDARD FEATURES**

- pendant control
- head end side rails, polypropylene head board & foot board
- battery back-up for emergency operation
- smooth, quiet operation
- fully welded and reinforced frame and deck
- beachwood powder coat finish
- width and length adjustable
- deck height adjustment from 16" low to 29.5" high
- model GEN1000 height adjustment 19.5" to 39"
- one year warranty all electrical components
- five year warranty on bed frame

Specifications subject to change without notice.  
Gendron Inc. Form No. PLBD4054(0109). All specifications and dimensions are approximate.