

GENDRON, INC.
 520 W Mulberry St
 Suite 100
 Bryan Ohio 43506

Effective January 1 2009

CUSTOMER SERVICE
 Phone: 800-537-2521
 Fax: 419-636-9261
 www.gendroninc.com

**EC+ Extra Care Bariatric Transport Chair
 Chair / Recliner / Stretcher**

**850 pound patient weight limit
 Model 6850**

Date:	P.O. #:
Dealer #:	Contact:
Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone No: (.....).....	Phone No: (.....).....

To order, check boxes on the price list below. Add up prices, insert "Order total." Mail, fax or call-in order.

BASE CHAIR MODEL NUMBER

6850 fully electric operation.....\$7600

OPTIONAL ACCESSORIES

Cushions

4", low shear cover, leg pad, seat & back ... Standard

Back Rest

electric operation Standard

Elevating Leg Rests

electric operation, independent operation .. Standard

Side Rails / Arm Rests

drop down, removable, chrome Standard

Push Handles

full width push handle Standard

IV Rod / Holder

322 stay lock IV holder\$210

Power Source

battery with on-board charger..... Standard

Caster Wheels

8" easy glide, four wheel central locking..... Standard

Convenience Items

W factory install weighing scales....consult factory

EC+ CHAIR AS SHOWN IN CATALOG

6850\$7600

STANDARD FEATURES

- 850 pound patient weight limit
- fully electric operation, pendant control
- one touch chair to stretcher conversion, 70" O.A length in stretcher position
- Four IV sockets
- 8" easy glide all surface casters
- Four wheel central lock and brake
- 4" high density cushions with low shear fabric cover
- full width push bar
- high / low 24.5" to 32.5" in stretcher position
- seat depth 22"
- seat width 29.5"
- back height 28.5"
- seat height 24.5", floor to top of seat cushion
- at widest point 33.5" wide
- approximate full cycles between battery charging: 12
- 40 – 50 second cycle for high / low articulation

Order total \$ _____

Specifications subject to change without notice.
 Gendron, Inc. PLST6850(0109)
 All specifications and dimensions are approximate.