

GENDRON[®]

TOTAL
SOLUTIONS
FOR
BARIATRIC
PATIENT
CARE

SAFETY.
COMFORT.
DIGNITY.

The Elite Acute Care Bed provides your bariatric patients with all of the above in an easy-to-use design, engineered to allow you to give your patients the best possible care.



Elite Bariatric Acute Care Bed




PATIENT HELPER/TRAPEZE

Transfer, lift and position patients with ease – swings 90 degrees from center for accessibility to the patient and eliminates the need to remove the trapeze from the bed when transferring the patient



CPR FOWLER BACKREST RELEASE

Respond faster to patient needs – lower backrest to flat within six seconds for CPR



FOUR WHEEL LOCK AND STEER

Steer safely and precisely – engage wheel lock and steer at any caster wheel



CONTROL PENDANT

Adjust settings confidently even in dark rooms – illuminated controls; five-function, self-storing at head end, left and right side



ADJUSTABLE DECK LENGTH

Ensure the comfort of your tall patients – one touch extension from 84" to 92" long



INTEGRATED BED SCALE SYSTEM

Ensure patient privacy – scale display window rotates to block the view of the information shown on the scale display



CAREGIVER CONTROL

Located at foot end of bed. Any or all hand pendant functions can be locked-out at the caregiver control to restrict patient access.

FEATURES INCLUDE:

- **BED POSITIONING/ ARTICULATION**
Provide proper clinical positioning from Trendelenburg and Reverse Trendelenburg to knee flex and cardio chair – easy-to-use button controls; one button cardio chair positioning and one button return to level positioning
- **OUT OF BED EXIT ALARM**
Alerts you when patient attempts to exit bed – programmable alarm sounds at the bed, with an option to connect to the nurse station




ADJUSTABLE DECK WIDTH

Transport safely and provide comfort – one touch extension from 39" to 48" to 54" wide



SIDE RAILS

Increase patient safety by preventing patient roll over – protect patients from falls with 17" high side rails



ELITE ACUTE CARE BED – MODEL 5154

FEATURES / SPECIFICATIONS



PATIENT WEIGHT CAPACITY: 1000 pounds

FOWLER BACK REST: 0-50 degrees

TRENDELENBURG/REVERSE TRENDELENBURG: 9 degrees

KNEE FLEX/GATCH: 0-20 degrees

DECK LOW POSITION: 17"

DECK HIGH POSITION: 26"

ADJUSTABLE DECK LENGTH: 84" - 92"

ADJUSTABLE DECK WIDTH: 39" - 48" - 54"

SIDE RAILS (4): 17" high above deck

CPR FOWLER BACKREST RELEASE: Head end, left and right

FOUR WHEEL LOCK AND STEER: Engage wheel locks or steer at any caster wheel

INTEGRATED PUSH HANDLES: Two, at head end

PATIENT HELPER/TRAPEZE: Swings 90 degrees from center

CONTROL PENDANT: Illuminated, five-function, self-storing

CAREGIVER CONTROL: Controls all functions and locks-out pendant functions

BED POSITIONING/ARTICULATION:

High/low, fowler backrest, Trendelenburg/Reverse Trendelenburg, knee flex/leg lift, cardio chair

ANGLE INDICATOR: Two at head end, one left, one right

UNDER BED ILLUMINATION: LED illumination

OUT OF BED EXIT ALARM: Sounds at bed with option to sound at nurse call station using DB 37 cable or 1/4" phone connector

INTEGRAL SCALE WEIGHING SYSTEM:

- Displays accurate weight in any articulated position including Trendelenburg and Reverse Trendelenburg and in any caster orientation
- One convenient control box contains both the scale controls and the caregiver key pad
- Scale display rotates – HIPAA-compliant
- Powered by bed electrical system and battery backup – does not require batteries and will function when disconnected from A/C power
- Weight change mode
- Dual weight display
- Freeze mode
- Backlit display
- Accurate to within 0.2% or 2 pounds, whichever is greater

WHY BUY GENDRON:

- 25 years experience serving the bariatric patient
- Designed and manufactured in the USA
- Complete bariatric product portfolio

GENDRON 

800.537.2521 • sales@gendroninc.com

www.gendroninc.com