


GENDRON®

Sentinel Bariatric Low Bed

SAFETY. COMFORT. DIGNITY.

The Sentinel Low Bed provides your bariatric patients with all of the above in an easy-to-use design, engineered to allow you to give your patients the best possible care.

ULTRA LOW DESIGN GIVES YOU AND YOUR PATIENTS PEACE OF MIND



STANDARD HEAD END SIDE RAILS
for patient safety

LOW DECK HEIGHT
Minimize potential for patient injury and caregiver risk due to accidental falls – lower deck to 9.25”

STURDY FRAME
for 850 lb. capacity



ADJUSTABLE SLEEP SURFACE WIDTH
Transport safely and provide comfort – one touch extension from 39” to 48” wide eliminates the need for multiple beds

MEETS ALL CRITICAL ARTICULATION
for clinical needs

TOTAL SOLUTIONS FOR BARIATRIC PATIENT CARE

SENTINEL BARIATRIC LOW BED – MODEL 3874 FEATURES / SPECIFICATIONS



PATIENT WEIGHT CAPACITY: 850 pounds

DECK LOW POSITION: 9.25"

DECK HIGH POSITION: 30"

DECK LENGTH: 86"

ADJUSTABLE DECK WIDTH: 39" - 48"

SIDE RAILS (2): 15" high above deck

CASTERS: Individual locking casters

CONTROLS: Hand pendant at head end

BATTERY POWER: Emergency battery is standard

CLINICAL ARTICULATION FUNCTIONALITY:

- High/low, fowler backrest
- Trendelenburg/Reverse Trendelenburg
- Knee flex/leg lift
- Cardio chair

WEIGHING SCALE: Available as factory installed option



WHY BUY GENDRON:

- 25 years experience serving the bariatric patient
- Designed and manufactured in the USA
- Complete bariatric product portfolio

GENDRON 

800.537.2521 • sales@gendroninc.com
www.gendroninc.com