

GENDRON, INC.

520 W Mulberry St
Suite 100
Bryan Ohio 43506

Effective January 1 2009

CUSTOMER SERVICE

Phone: 800-537-2521
Fax: 419-636-9261
www.gendroninc.com

**Series 1190 Bariatric Extra Heavy Duty Stretcher 1000 pound patient weight limit
Transport, Transfer**

Date: P.O. #:

Dealer #: Contact:

Bill To: Ship To:

Address: Address:

City/State/Zip: City/State/Zip:

Phone No: (.....)..... Phone No: (.....).....

To order, check boxes on the price list below. Add up prices, insert "Order total." Mail, fax or call-in order.

BASE STRETCHER MODEL NUMBER

1190 Powder coat* frame and top\$4430

OPTIONAL ACCESSORIES

Litter Pads

4" hinged litter pad (high density) Standard
 324-13 replacement litter pad.....\$425

Fowler Back Rest Attachment

132 Hydraulic, powder coat* finish Standard

Side Rails

Three position vertical, chrome Standard

Utility Shelf

140 Full length, powder coat* finish.....\$150

IV Rod / Holder

121 Chrome.....\$55
 122 Chrome, telescopic.....\$75
 136 Stainless steel, telescopic\$80
 139 Chrome telescopic for pump mounting\$100

Oxygen Holder

144 Powder coat* finish.....\$80

Caster Wheels

165 Central lock and brake..... Standard

Convenience Items

119 Restraint straps\$85 ea.
 130 Padded arm board.....\$175
 172 Chart holder.....\$125

STRETCHER AS SHOWN IN CATALOG

1190\$4430
Pictured with 4" hinged pad, Fowler back rest, three position vertical side rails

STANDARD FEATURES

- Welded steel frame, 1000 pound patient weight limit
- Litter top, 80" L x 35" W
- Three position, vertical side rails
- Four IV sockets
- Full perimeter non-marring bumper
- 8" all surface casters
- Four wheel central lock and brake
- Beechwood powder coat* epoxy finish

Order total \$ _____

Specifications subject to change without notice.
Gendron, Inc. PLST1190(0109)
All specifications and dimensions are approximate.