

GENDRON, INC.

520 W Mulberry St
Suite 100
Bryan Ohio 43506

Effective January 1 2009

CUSTOMER SERVICE

Phone: 800-537-2521
Fax: 419-636-9261
www.gendroninc.com

**EC+ Extra Care Bariatric Stretcher
ER, Trauma, Transport**

**1000 pound patient weight limit
Model 1600**

Date:
Dealer #:
Bill To:
Address:
City/State/Zip:
Phone No: (.....).....

P.O. #:
Contact:
Ship To:
Address:
City/State/Zip:
Phone No: (.....).....

To order, check boxes on the price list below. Add up prices, insert "Order total." Mail, fax or call-in order.

BASE STRETCHER MODEL NUMBER

1600 fully electric operation.....\$5600

OPTIONAL ACCESSORIES

Litter Pads

4" hinged litter pad, low shear cover..... Standard

Fowler Back Rest Attachment

electric operation Standard

Side Rails

full length fold down, chrome..... Standard

Push Bar

push bar at head end Standard

IV Rod / Holder

322 stay lock IV holder\$210

Power Source

battery with on-board charger..... Standard

Caster Wheels

8" easy glide, four wheel central locking..... Standard

Convenience Items

W factory install weighing scales....consult factory
 119 Restraint straps\$85ea

STRETCHER AS SHOWN IN CATALOG

1600\$5600

STANDARD FEATURES

- 1000 pound patient weight limit
- fully electric operation, pendant control
- fold down side rails
- Four IV sockets
- four corner non-marring bumpers
- 8" easy glide all surface casters
- Four wheel central lock and brake
- 4" high density litter pad with low shear fabric cover
- integral self storing push handles
- CPR release, head end, left & right side access
- high / low 20" to 28"
- fowler angle 55 degrees
- knee gatch angle 48 degrees
- trendelenburg / reverse trendelenburg 10 degrees
- patient surface 37" x 80"
- overall 39" x 83"
- approximate full cycles between battery charging: 12

Order total \$ _____

Specifications subject to change without notice.
Gendron, Inc. PLST1600(0109)
All specifications and dimensions are approximate.