

GENDRON, INC.
 520 W Mulberry St
 Suite 100
 Bryan Ohio 43506

CUSTOMER SERVICE
 Phone: 800 - 537 - 2521
 Fax: 419-636-9261
 www.gendroninc.com

Bariatric Seating
Weight Limit: 750 to 850 Pounds

Price List Effective January 1 2009

Date
 Dealer Name:
 Account Number:
 Address:
 City/State/Zip:
 Phone No.: (.....).....
 Authorized Signature

P.O. #:
 Contact:
 Ship To:
 Address:
 City/State/Zip:
 Phone No.: (.....).....
 Therapist / Referral:

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, fax or call-in order.

MODEL NUMBER

- 7125 25"W x 20"D x 18" seat height \$1045
- 7130 30"W x 20"D x 18" seat height \$1100
- 7140 40"W x 20"D x 18" seat height \$1135
- 7126 25"W x 20"D x 18" – 23" adjustable seat height.... \$1175
- 7131 30"W x 20"D x 18" – 23" adjustable seat height ... \$1225
- 7141 40"W x 20"D x 18" – 23" adjustable seat height.... \$1260

Arm Style

- full length..... standard

Back Style

- permanent.. standard

Features include:

- seat height and sit to stand position available with adjustable seat height
- arms & legs are laminated hardwood
- front and rear rails are mortise & tenon joined and glued
- seat and back platform are removable
- arms and legs are natural wood finish
- seat & back upholstery color Marina blue
- seat and back fabric meets CAL 133 code
- seat width measured between arms
- add 4" to seat width to determine overall width
- 750 pound patient weight capacity

BARIATRIC PATIENT ROOM RECLINER

- 7150 patient recliner 850 pound capacity\$2395
- 7155 patient recliner 850 pound capacity\$2185

Features include:

- attendant operated recline to 40 degrees
- leg rest articulates when back is lowered / reclined
- upholstered high density foam for user comfort
- model 7150 seat width 29.25" / model 7155 seat width 25.25"
- seat depth 19"
- navy, portwine, brown, grey color choices
- seat to floor 20.25"
- easily portable, swivel and locking casters
- model 7150 overall maximum width 38.25"
- model 7155 overall maximum width 34.25"

All specifications and dimensions are approximate.
 Specifications subject to change without notice.
 Gendron Inc. Form No. PLSTSeating (0109).

