

GENDRON, INC.
520 W Mulberry St.
Suite 100
Bryan Ohio 43506

Effective January 1 2009

CUSTOMER SERVICE
Phone: 800-537-2521
Fax: 419-636-9261
www.gendroninc.com

Maxi Rest Bariatric Beds

Weight Capacity: 750 Pounds

Date: P.O. #:
Dealer #: Contact:
Bill To: Ship To:
Address: Address:
City/State/Zip: City/State/Zip:
Phone No: (.....) Phone No: (.....)

MODEL NUMBER

- 4748B 48 x 80 sleep surface, 3 function electric, one piece frame, head end rails, head & foot boards..... \$4995
- 4748D 39 x 80 sleep surface, 3 function electric, one piece frame, head end rails, head & foot boards, width and length adj. \$5925
- 4748DX 39 x 80 sleep surface, 5 function electric, one piece frame, head end rails, head & foot boards, width and length adj. \$7410
- 4748SB 48 x 80 sleep surface, 3 function electric, two piece frame, head end rails, head & foot boards..... \$5095
- 4748SD 39 x 80 sleep surface, 3 function electric, two piece frame, head end rails, head & foot boards, width and length adj. \$6240
- 4748SDX 39 x 80 sleep surface, 5 function electric, two piece frame, head end rails, head & foot boards, width and length adj. \$7775

Note: three function electric features height adjustment, fowler back rest and leg lift
five function electric features height adjustment, fowler back rest, leg lift, trendelenburg and reverse trendelenburg

Important: models 4748D, 4748SD, 4748DX and 4748SDX adjust from 39" to 48" wide utilizing integral extension panels, length adjustable to 88"

OPTIONAL ACCESSORIES

Mattress

- 3980F 6" bariatric mattress, 39 x 80 \$820
- 4880F 6" bariatric mattress, 48 x 80 \$890
- 4880-1F 6" bariatric mattress, 39-48 x 80 \$1150
- 5886F 6" bariatric mattress, 48 x 84 850 lb. cap... \$1040

Note: model 4880-1F mattress includes side bolsters to expand to 48" wide

Bolsters / Foot End

- 3908F 6"x8"x39" mattress bolster with deck panel..... \$350
- 4808F 6"x8"x48" mattress bolster with deck panel \$375
- 22155 deck extension panel only \$250

Side Rails

- head end side rails, pair, included with bed N/C.
- 4043 one pair side rails, foot end for D / DX series ... \$565
- 4043-1 one pair side rails, foot end for "B" series bed. \$450
- 4044-1 one pair side rails extra tall for D / DX beds.... \$215

Trapeze Frame

- 4045 over bed trapeze and frame 80" long \$1250
- 4046 over bed trapeze and frame 88" long \$1450

Accessory Wall Bumper

- 4047 head or foot end full width wall bumper \$215

Scale

- W integral weighing scales, caster stem mount . \$5000

Caster Wheels

- 5" hospital casters, two locking Standard

CPR Release

- CPR fowler back rest release \$375

NOTE: CPR release is not available on split frame beds

STANDARD FEATURES

- pendant control
- head end side rails, polypropylene head board & foot board
- battery back-up for emergency operation
- smooth, quiet operation
- fully welded and reinforced frame and deck
- beachwood powder coat finish
- width and length adjustment expands to 48" wide and 88" long
- Maxi Rest deck height adjustment from 15.5" low to 29.5" high
- one year warranty all electrical components
- five year warranty on bed frame

Specifications subject to change without notice.
Gendron Inc. Form No. PLBD4748(0109)
All specifications and dimensions are approximate.