

GENDRON, INC.  
 520 W Mulberry St  
 Suite 100  
 Bryan Ohio 43506

CUSTOMER SERVICE  
 Phone: 800.537.2521  
 Fax: 419-636-9261  
 www.gendroninc.com

Price List Effective January 1 2009

**Regency XL 2002 Baritaric Manual Wheelchairs  
 Seat Widths: 24" 26" 28"**

**Models 5624 5626 5628  
 Weight Limit: 600 Pounds**

Date .....  
 Dealer Name: .....  
 Account Number: .....  
 Address: .....  
 City/State/Zip: .....  
 Phone No.: (.....) .....  
 Authorized Signature .....

P.O. #: .....  
 Contact: .....  
 Ship To: .....  
 Address: .....  
 City/State/Zip: .....  
 Phone No.: (.....) .....  
 Therapist / Referral: .....

To order, check boxes on the price list below. Add up prices, insert "total chair price", and mail, [e-mail](#), fax or call-in order.

**BASE WHEELCHAIR MODEL NUMBER**

56 Series base price ..... \$ 1255

**Seat Dimensions W x D inches**

5624-18 24 x 18 ..... N/C  
 5626-18 26 x 18 ..... N/C  
 5628-18 28 x 18 ..... N/C

**Seat to Floor (STF) Height**

(9) 19.5" ..... N/C

**Arm Style**

(-61) desk length, removable, 8" fixed height ..... standard  
 (-62) full length, removable, 8" fixed height ..... N/C  
 (-63) full length, removable, 8"-12" height adjustable ..... \$150

**Back Style**

fixed height 16" with 5 degree back angle ..... standard

**Frame Colors**

Chrome powder coat ..... standard

**Front Riggings**

(-33) swing away footrests ..... \$120  
 (-36) swing away elevating legrests ..... \$235

**Wheel Options**

24" mag composite wheels ..... standard

**Caster Options**

8" x 1-1/4" semi pneumatic wheels ..... standard

**Seat & Back Cushions**

75 padded solid back insert ..... \$160  
 94 padded solid seat insert ..... \$160

**IV & Oxygen Tank Holder**

(-12) oxygen holder ..... \$135  
 (-95) IV holder/height adjustable rod ..... \$105  
 (-1295) oxygen tank/IV holder combination ..... \$170

**Miscellaneous Accessories**

(-59) leg guard ..... \$100  
 (-80) safety belt, specify  auto buckle  Velcro ..... \$80  
 (-83) foot plate heel loops ..... \$60  
 (-85) anti tippers, rear with wheels ..... \$125  
 (-86) wheel lock extensions ..... \$50  
 (-93) patient records holder ..... \$55  
 (-STNCL) custom name stenciling ..... \$75

**Order Total \$ \_\_\_\_\_**

**Standard features include:**

- heavy duty steel frame, 5/8" rear axle
- heavy duty cross-brace
- cushioned arm rests with steel clothing guards
- through bolted seat upholstery
- black nylon upholstery with heavy duty seat liner
- gussets reinforced & welded caster fork & caster journal
- back angle 10 degrees
- 8" x 1-1/4" semi-pneumatic caster wheels
- 24" mag composite wheels with molded hand rims

All specifications and dimensions are approximate  
 Specifications subject to change without notice.  
 Gendron Inc. Form No.: PLWCXL2002 Series (0109)