

GENDRON, INC.
 520 W Mulberry St
 Suite 100
 Bryan Ohio 43506

Effective January 1 2009

CUSTOMER SERVICE
 Phone: 800-537-2521
 Fax: 419-636-9261
 www.gendroninc.com

Specialty Beds
Extra Care Bariatric Low Bed / Extra Care Ultra Low Bed

Date: P.O. #:
 Dealer #: Contact:
 Bill To: Ship To:
 Address: Address:
 City/State/Zip: City/State/Zip:
 Phone No: (.....) Phone No: (.....)

To order, check boxes on the price list below. Mail, fax or call-in order.

MODEL NUMBER

- 3607 36 x 80 sleep surface, 3 function electric, head end rails, head & foot boards, width and length adj. \$6225
- 3848 48 x 84 sleep surface, 5 function electric, head end rails, head & foot boards, height 9" - 26.75" \$6700

Note: three function electric features height adjustment, fowler back rest and leg lift
 five function electric features height adjustment, fowler back rest, leg lift, trendelenburg and reverse trendelenburg

Important: model 3607 features expandable bed deck width, refer to optional width extension kits shown below
 model 3607 height range: low position 7" to high position 28"

Patient Weight Capacity: Model 3607 / 500 pounds Model 3848 / 850 pounds

OPTIONAL ACCESSORIES

Mattress

- 3608F** 6" bariatric mattress, 36 x 80 \$820
- 4280F** 6" bariatric mattress, 42 x 80 \$890
- 5886F** 6" bariatric mattress, 48 x 84 850 lb. cap. . \$1040

Mattress Bolsters

- 3908F** 4"x8"x42" mattress bolster foot end \$350
- 4808F** 4"x8"x48" mattress bolster foot end \$350
- 0680F** 6"x6"x80" mattress bolster side pair \$425
- 3604-1F** 4"x6"x36" foot end bolster \$200

Side Rails

- head end side rails, pair, included with bed N/C.
- 4043 one pair side rails, foot end model 3848 \$565
- 4043-4 one pair side rails, foot end model 3607 \$950

Deck Extension Kits for Ultra Low

- 14203 kit to expand deck width to 42" \$500
- 14806 kit to expand deck width to 48" \$590

CPR Release

- CPR fowler back rest release \$375

Accessories

- 13604 foot end deck extender 4"x36" for model 3607 . \$225
- 10121 IV rod / holder chrome \$100
- Attendant push bar models 3607 or 3848 Standard

STANDARD FEATURES

- pendant control
- head end side rails, polypropylene head board & foot board
- battery back-up for emergency operation
- smooth, quiet operation
- fully welded and reinforced frame and deck
- beachwood powder coat finish
- one piece bed frame models 3607 and 3848
- one year warranty all electrical components
- five year warranty on bed frame

Specifications subject to change without notice.
 Gendron Inc. Form No. PLBDSPEC(0109)
 All specifications and dimensions are approximate.

Order total \$ _____