

**GENDRON, INC.**  
 520 W Mulberry St.  
 Suite 100  
 Bryan Ohio 43506

**CUSTOMER SERVICE**

Phone: 800-537-2521  
 Fax: 419-636-9261  
 Email: [www.gendroninc.com](http://www.gendroninc.com)

**EFFECTIVE DATE:** March 2018

**NOTE:** Specifications subject to change without notice.  
 Gendron Inc. Form No. PLWC 6800 Series(0318)  
 All specifications and dimensions are approximate.

**PRODUCT:** Regency XL 2000 Bariatric Manual Wheelchairs  
 Model Series 6800  
 Seat Widths: 26" 28" 30" 32" 34"  
 Weight Limit: 850 Pounds

**BASE WHEELCHAIR MODEL NUMBER**      **DESCRIPTION**  
 68      Base wheelchair

**Seat Dimensions W x D Inches**

6826-20	26 x 20
6828-20	28 x 20
6830-20	30 x 20
6832-20	32 x 20
6834-20	34 x 20
6826-22	26 x 22
6828-22	28 x 22
6830-22	30 x 22
6832-22	32 x 22
6834-22	34 x 22

**Seat to Floor Height**

(8)	18" (1) Adult
(7)	17" (1) Adult
(6)	16" (1) Hemi
(5)	15" (2) Super hemi (includes 20" and 6" wheels)
(1)	24" Rear / 8" caster standard for 16", 17", 18" seat hts.
(2)	20" Rear / 6" caster required for 15" seat ht.

**Arm Style**

(-62)	Full length, removable
(-63)	Full length, removable 8" - 12" height adjustable

**Push Bar**

Full width removable push bar

**Back Style**

Adjustable height 16", 17", 18" (factory set at 16")  
 Power assist 45 degree reclining back, 36" high

**Note:** Push bar & rear anti tips are standard on recliners

**Frame Color**

Black

**Front Riggings**

(-26)	Adjustable front footrests
(-27)	Adjustable front elevating legrests

**DESCRIPTION**

**Wheel Options** 24" Metal wheels solid tires for 16", 17", 18" height  
 20" Mag wheels with semi pneumatic 15" height  
**Note:** Wheels have steel hand rims

**Caster Options** 8" x 2" Semi pneumatic wheels for 16", 17", 18" height  
 6" x 2" Semi pneumatic wheels for 15" height

**IV & Oxygen Tank Holder**

(-12)	Oxygen holder
(-13)	Oxygen holder for recliner
(-14)	IV rod/holder for recliner
(-95)	IV holder height adjustable
(-1295)	Oxygen tank/IV holder combination

**OPTIONAL ACCESSORIES**

(-80)	Safety belt, specify __Auto buckle __Velcro
(-85)	Anti tippers, rear with wheels

**STANDARD FEATURES INCLUDE**

Heavy duty steel frame, 5/8" rear axle  
 Double cross-braces  
 Cushioned arm rests with steel clothing guards  
 Through bolted seat upholstery  
 Black nylon upholstery with heavy duty seat liner  
 Caster fork & caster journal-gussets reinforced & welded  
 Back angle 10 degrees  
 8" x 2" Semi pneumatic caster wheels  
 24" Metal mag wheels with steel hand rims  
 Adjustable seat to floor height  
 All models feature heavy-duty adjustable front riggings  
 Full length push bar across back