

**GENDRON, INC.**

520 W Mulberry St.  
Suite 100  
Bryan Ohio 43506

**CUSTOMER SERVICE**

Phone: 800-537-2521  
Fax: 419-636-9261  
Email: [www.gendroninc.com](http://www.gendroninc.com)

**EFFECTIVE DATE:** March 2108

**NOTE:** Specifications subject to change without notice.  
Gendron Inc. Form No. PLWC67(0318)  
All specifications and dimensions are approximate.

**PRODUCT:** Regency XL 2000 Bariatric Manual Wheelchairs  
Model Series 6700  
Weight Limit 700-750 Pounds  
Seat Widths: 18" 20" 22" 24" 26" 28" 30" 32" 34"

<u>BASE MODEL NUMBER</u>	<u>DESCRIPTION</u>
67	Base wheelchair frame
<b>Seat Dimensions W x D Inches</b>	
6730-18	30 x 18
6732-18	32 x 18
6734-18	34 x 18
6720-20	20 x 20
6722-20	22 x 20
6724-20	24 x 20
6726-20	26 x 20
6728-20	28 x 20
6730-20	30 x 20
6732-20	32 x 20
6734-20	34 x 20
6720-22	20 x 22
6722-22	22 x 22
6724-22	24 x 22
6726-22	26 x 22
6728-22	28 x 22
6730-22	30 x 22
6734-22	34 x 22
(-18)	18" (1) Adult 700 pound capacity
(-17)	17" (1) Adult 700 pound capacity
(-16)	16" (1) Hemi 700 pound capacity
(-15)	15" (2) Super hemi (includes 6" & 20" wheels 700 lb)
(1)	24" Rear / 8" caster standard for 16", 17", 18" seat hts.
(2)	20" Rear / 6" caster required for 15" seat ht.
	Wheel Package Up-Grade for 750 pound weight capacity
	24" Metal mag wheels pair, for 16", 17", 18" seat hts.
	20" Mag wheels pr./ 6" x 2" caster wheels pr. for 15" seat ht.
<b>Arm Style</b>	
(-61)	Desk length, removable, 8" fixed ht.
(-62)	Full length, removable, 8" fixed ht.
(-63)	Full length, removable, 8"-12" ht. adjustable
<b>Back Style</b>	
(-78)	Adjustable ht. 16", 17", 18" (factory set at 16") Power assist 45 degree reclining back, 36" high
	<b>Note:</b> Push bar and rear anti tips are standard on recliners
<b>Frame Colors</b>	
	Black
<b>Front Riggings</b>	
33-1	S/A footrests up to 600 lbs cap
36-1	S/A elevating legrests up to 600 lbs cap
28-1	S/A foot rests over 600 lbs cap
29-1	S/A elevating legrests over 600 lbs cap

<u>Wheel Options</u>	<u>DESCRIPTION</u>	<u>BASE PRICE</u>
	24" Uni mold mag wheels, pair, for 16", 17", 18" seat hts.	
	20" Mag wheels semi pneumatic pair, for 15" seat ht.	
<b>Caster Options</b>		
	8" X 2" Semi pneumatic wheels for 16", 17, 18" seat hts.	
	6" X 2" Semi pneumatic wheels for 15" seat ht.	
<b>IV &amp; Oxygen Tank Holder</b>		
(-12)	Oxygen holder	
(-13)	Oxygen holder for recliner wheelchair	
(-14)	IV holder for recliner wheelchair	
(-95)	IV holder/ht. adjustable rod	
(-1295)	Oxygen tank/IV holder combination	

<u>OPTIONAL ACCESSORIES</u>	
(-80)	Safety belt, specify __Auto buckle __Velcro
(-85)	Anti tippers, rear with wheels
(-40)	Full width push bar