

**GENDRON, INC.**  
 520 W Mulberry St.  
 Suite 100  
 Bryan Ohio 43506

**EFFECTIVE DATE:** March 2018

**PRODUCT:** Regency XLC Wheelchairs Custom Bariatric Manual  
 and Power Wheelchairs Seat Widths: 20" Plus  
 Weight Limit: Up to 1000 Pounds

**RIGID FRAME WHEELCHAIR BASE**

7850 Manual base  
 7800 Power base  
 Power base right control \_\_\_\_ Power base left control \_\_\_\_

Seat Width		
20"	Standard	31"
21"	\$425	32"
22"	\$575	33"
23"	\$700	34"
24"	\$825	35"
25"	\$1,025	36"
26"	\$1225	37"
27"	\$1425	38"
28"	\$1600	39"
29"	\$1775	40"
30"	\$2025	*Call for seat widths greater than 40"

**Note:** During manufacture - rear wheels on power base chairs can be "tucked" under seat frame.  
 Consult factory for narrowest overall width based on seat dimension.

**Seat Construction**  
 Solid pan, no cushion, metal deck  
 Solid with 2" foam  
 Solid with 4" foam

**Back Style**  
 Permanent  
 Manual recliner 45 degrees \_\_\_\_ 60 degrees \_\_\_\_

**Power Accessories / Back / Tilt in Space**  
 Power recliner 45 degrees \_\_\_\_ 60 degrees \_\_\_\_  
 Power recliner 90 degrees

**Back Support Surface**  
 Solid back, no padding  
 Solid with 2" foam pad  
 Solid with 4" foam pad

**Arm Style**  
 Full length, height adjustable, 7" - 9.5", removable  
 Full length, height adjustable, 9" - 11", removable

**Front Riggings/Heavy Duty/Reinforced**  
 Footrest for seat height 14" - 16"  
 Footrest for seat height 16" - 20"  
 Elevating legrests for seat height 14" - 16"  
 Elevating legrests for seat height 16" - 20"

**CUSTOMER SERVICE**

Phone: 800-537-2521  
 Fax: 419-636-9261  
 Email: [www.gendroninc.com](http://www.gendroninc.com)

**NOTE:** Specifications subject to change without notice.  
 Gendron Inc. Form No. PLWC RegencyXLC(0318)  
 All specifications and dimensions are approximate.

**Frame Colors**  
 Black Standard

Seat Depth		
18"	Standard	23"
19"	\$575	24"
20"	\$625	25"
21"	\$675	26"
22"	\$850	

**Seat Height**  
 14" & 15" available on manual wheelchairs only  
 16", 17", 18", 19" manual or power wheelchairs  
 20" Manual or power wheelchairs  
 21" Manual or power wheelchairs  
 \*Measured from floor to seat pan without foam padding.  
 For any other seat heights/depths not listed, call customer service.

**Seat Plane / Angle**  
 3 Degrees  
 Other seat angles, specify distance from floor for:  
 Front edge of seat \_\_\_\_ Rear edge of seat \_\_\_\_

**Back Height**  
 16"  
 17" \_\_\_\_ 18" \_\_\_\_ 19" \_\_\_\_ 20" \_\_\_\_  
 21" \_\_\_\_ 22" \_\_\_\_ 23" \_\_\_\_ 24" \_\_\_\_  
 25" \_\_\_\_ 26" \_\_\_\_ 27" \_\_\_\_ 28" \_\_\_\_

**Back Angle**  
 Back angle 7 degrees  
 Other back angle, specify \_\_\_\_

**OPTIONAL ACCESSORIES**

- (-12) Oxygen tank holder (see Note)
- (-80) Safety belt, specify \_\_\_\_Auto buckle \_\_\_\_Velcro
- (-85) Wheeled anti-tip brackets / rear
- (-86) Wheel lock extensions, manual chairs
- (-95) IV rod and holder (see Note)
- (-1295) Oxygen tank/IV rod holder combination (see Note)
- P20529 Battery 12V, gel cell

**Note:** IV/oxygen holder not available on recliners & tilt in space