

MAXI CARE BARIATRIC LIFT

The Gendron Maxi Care Bariatric Lift was designed and built with the safety of the patient and care giver in mind.

The Maxi Care Bariatric Lift is fully electric and the swivel and lower functions dramatically reduce the risk of caregiver injury when transferring a patient. The stability of the base legs open and extend to ensure a safe patient transfer.



SPECIFICATIONS

Capacity:	2 Models - 700# (317.5 kg) and 1000# (453.6 kg)
Configuration of Lifting Mechanism:	Battery powered linear actuators
Boom - Pivot Height x Overall Length:	58" x 63" (147.3cm x 1600.2cm)
Patient Support:	Various Slings with 4-point suspension
Lift Range (lowest position):	Floor
Lift Range (highest position):	20½" (52.1cm) to bottom of sling
Base Size L x W - transport & storage (Legs retracted & closed):	53" x 36" (134.6cm x 58.4cm)
Base Size L x W - in operation (Legs extended & open):	58" x 82½" (147.3cm x 209.6cm)
Note: Leg movement is electronically controlled	
Minimum Under-bed Clearance Requirement:	4½" (15.2cm)
Minimum Ceiling Height Requirement:	86" (218.4 cm)
Brakes:	Rear Locking Casters
Battery Type:	Rechargeable (Integrated Charging System)/ Removable/Replaceable
Operating Controls:	Caregiver operated, hand-held, pushbutton pendant with Manual lower capability

GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
Phone: 419.636.0848 www.gendroninc.com