



MAXI REST BARIATRIC BED

MAXI REST 4842 SERIES

Features

- Robust steel two-piece split-frame construction for ease of delivery and transport
- Bed deck width is expandable from 39" to 42" to 48"
- Bed deck length is adjustable from 80" to 88"
- 5" Deluxe Caster Set includes 2 Total-Lock and 2 Steer-Lock for ease of maneuverability
- Standard 3-Function Positioning hand control
- Optional 5-Function Positioning hand control (Trendelenburg / Reverse Trendelenburg and Chair positioning)
- Synchronized Hi/Lo Motors ensure bed remains level
- Battery backup standard for dependability
- Optional Intuitive side rail system for both head and foot sections
- Optional laminate wood grain or plastic molded head and foot board options
- Optional factory-installed weight scale system
- Safe working load: 800 lb (362.9 kg)*



* Refer to user manual

Dimensions & Specifications

| Bed Model | Patient Weight Capacity EVENLY DISTRIBUTED | | Height Range | | overall Length | | Sleep Deck Width Range | | Number of Functions | Rolling Position | Number of Motors | Built-in Scale | Emergency Battery Backup |
|-----------|---|-------|--------------|-----------|----------------|----------------|------------------------|---------------------|-------------------------|------------------|------------------|----------------|--------------------------|
| | lb | kg | inch | cm | inch | cm | inch | cm | | | | | |
| 4842 | 800 | 362.9 | 13.5-27.5 | 34.3-69.9 | 80 or 88 | 203.2 or 223.5 | 39-42 or 48 | 99.1-106.7 or 121.9 | 3 (std) or 5 (optional) | height | four | optional | battery |

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



MAXI REST 4842 SERIES

Maxi Rest 4842 Accessories

- GL928000.....38" Laminate (set of head and foot boards)
- GL929000.....39" HDPE Plastic (set of head and foot boards)
- GL930000.....48" HDPE Plastic (set of head and foot boards)
- GL927000.....Standard Chrome Rails (pair)
- GL927020.....Standard Chrome Rails w/ Pull-outs for 39", 42", 48"
- GL93100.....Trapeze 80" Frame Mount Full Length
- GL93180.....Trapeze 88" Frame Mount Full Length
- GL93200.....Trend/Rev Trend Hand Control (5-function)
- GL93300.....Wall Bumper Full Width
- GL93600L.....Egress Bar (Assist Bar) and Left Hand Mounting Pullout (48" width only)
- GL93600R.....Egress Bar (Assist Bar) and Right Hand Mounting Pullout (48" width only)
- GL94000.....8" Bolster Mattress and Ext. Pan Kit

Warranty


- Main Frame and welds:..... Five years
- Control Box and Actuator Motors:..... Two years
- Hand Control Pendant, Staff Control, and Cabling:..... Two years
- Head Board and Foot Board:..... One year
- All other durable components not listed above:..... Two years




39" HDPE Plastic (set of head and foot boards)



Standard Chrome Rails (pair),
shown with optional 5-function Hand Control



Trapeze 80" Frame Mount Full Length,
Deluxe Plastic Rails (set- LH/RH)



Trend/Rev Trend Hand Control
(5-function)



Optional factory-installed scale system

WARNING: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com