



WHEELCHAIRS

MRI TRANSPORT WHEELCHAIR

The Gendron MRI Transport Wheelchairs are designed specifically for use in the MRI suite, making it ideal for hospitals and medical centers everywhere. These transport wheelchairs are constructed of non-magnetic materials making it safe for use in the suite, however they should not be utilized directly inside of the MR system. Proudly made in the USA, the Gendron MRI Transport Wheelchair is designed with a distinctive blue color (4650 and 4850 models only) to help staff immediately identify MR safe wheelchair equipment. Each chair features swing away footrests as well as removable full length armrests for comfort and convenience. Ideal for transferring patients from one department to another, these wheelchairs are available in 350 lb, 650 lb and 850 lb weight capacities.

Special Notice:

Non-clinical testing demonstrated that each product (4650MR / 4850MR) is MR conditional and can be used in the MRI environment according to the following conditions:

Static magnetic field of 3 - Tesla or less

Spatial gradient magnetic field of 720 - Gauss/cm or less

Note: these products are intended for use inside the MRI environment (e.g., in the MR room, close to the scanner). However, they should not be utilized directly inside of the MR system (e.g., inside the bore of the scanner).

Features

- Made with materials that are non-magnetic for use inside the MR suite
- Available in 20" seat width, 18" seat depth and low 18" seat height for 350 lb weight capacity model
- Available in 26" seat width, 20" seat depth and low 18" seat height for both 650 and 850 lb weight capacity models
- Black, nylon upholstery
- Removable full length arm style with padded armrests
- Adjustable standard swing away footrest with durable black footplates
- Rear 24" x 1" wheels with black handrims
- Rear anti-tippers included
- Adjustable back height from 16" or 18"
- Maximum Weight Capacity 350 lb - EVENLY DISTRIBUTED on Model 4000MRQ2
- Maximum Weight Capacity 650 lb - EVENLY DISTRIBUTED on Model 4650MR and 850 lb on 4850MR



4650/4850 model shown

Warranty

One (1) year limited warranty

Specifications:

	Seat Width		Seat Depth		Seat Height		Back Height from seat		Rear Wheel diameter		Front wheel diameter		overall Height		overall Width		overall folded Width		shipping Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg			
4000MRQ2	20	50.8	18	45.72	18	45.72	16	40.64	24	60.96	8	20.32	36	91.44	35	88.9	21	53.34	81	36.74	350	158.76	1/ea
4650MR	26	40.64	20	50.8	18	45.72	16	40.64	24	60.96	8	20.32	36	91.44	35	88.9	21	53.34	81	36.74	650	294.84	1/ea
4850MR	26	40.64	20	50.8	18	45.72	16	40.64	24	60.96	8	20.32	36	91.44	35	88.9	21	53.34	81	36.74	850	385.55	1/ea

WARNINGS: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an ISO 13485:2016 MDSAP Certified Company

© 2020, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com