



SHOWER CHAIR / COMMODE

BARIATRIC COMMODE

Features

- 750 lb maximum weight capacity (Evenly Distributed)
- Powder-coated steel frame
- Adjustable, removable armrest for lateral transfer
- Multiple Seat width options
- Pail is removable from the back



Specifications

Bariatric Commode Model	Patient Weight Capacity EVENLY DISTRIBUTED		Overall Width		Seat Width		Seat Depth		Seat to Floor Height		Frame Material	Arm Style
	lb	kg	inch	cm	inch	cm	inch	cm	inch	cm		
523326	750	340.2	31.4	79.8	26	66.0	20	50.8	20 - 24	50.8 - 61.0	Coated Steel	Removable / Adjustable
523330	750	340.2	31.4	79.8	30	76.2	20	50.8	20 - 24	50.8 - 61.0	Coated Steel	Removable / Adjustable

WARNING: Cancer and Reproductive Harm -
www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

