



BARIATRIC HOMECARE BED



3648 SERIES

**Features**

- Two-piece frame construction for ease of delivery and transport
- Solid pan-style sleep deck
- Deluxe caster set features durable Total-Lock and Steer-Lock casters
- Fully electric, pendant control
- Standard 3-Function Positioning hand control
- Maximum weight capacity, EVENLY DISTRIBUTED: 650 lb (194.8 kg)
- Head end side rails (shown on bed at right) are included



**Dimensions & Specifications**

Bed Model	Patient Weight Capacity EVENLY DISTRIBUTED		overall Length		Sleep Deck Width	
	lb	kg	inch	cm	inch	cm
3648	<b>650</b>	<b>194.8</b>	84	213.4	48	121.9



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2021, GF Health Products, Inc. All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.