



TRAPEZE / FLOOR STANDS

PORTABLE PATIENT HELPER TRAPEZE

The Gendron® Portable Patient Helper Trapeze unit is an important patient room accessory designed to help the users change positions while in bed and to aid in transferring from bed to chair with easy attendant assistance.

Features

- Durable, robust steel frame
- Helps to position patient in bed and transfer from bed to other equipment with minimal assistance
- Handle bar location fully adjustable using strap with buckle to maximize patient accessibility
- Comes with 2" diameter x 1" wide casters and floor stand base that folds for transportability
- Limited 2 Year Limited Warranty



Specifications:

| | overall Height | | overall Width (frame and leg width) | | overall Depth (leg length) | | Boom Length (extends horizontally over bed) | | Product Weight | | Weight Capacity EVENLY DISTRIBUTED | | UOM | |
|--|----------------|----|--|----|-------------------------------|----|--|----|----------------|-----|---------------------------------------|-----|--------|------|
| | inch | cm | inch | cm | inch | cm | inch | cm | lb | Kg | lb | Kg | | |
| | 5160 | 72 | 182.88 | 36 | 91.44 | 63 | 160.02 | 35 | 88.90 | 154 | 69.85 | 850 | 385.55 | 1/ea |

WARNINGS: Cancer and Reproductive Harm - www.P66Warnings.ca.gov

Color variations may occur due to monitor or printer settings.

GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
MEDICAL DEVICE SINGLE AUDIT PROGRAM
Intertek

+1.770.368.4700
Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

graham field

GF200026RevA20

© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.