



MATTRESSES

The Bariatric Care Foam Mattresses offer comfort, support and durability for the bariatric patient.

Features:

- Super soft tapered heel section provides added pressure reduction for sensitive heel and ankle
- Soft and durable top covers are bacteria resistant, flame retardant, fluid resistant, anti-fungal, anti-static, non-allergenic, have near 0% vapor transmission, and are easy to clean
- Soft and durable Vinyl bottom covers are bacteria resistant, flame retardant, fluid resistant, anti-static, non-allergenic, stain resistant, and easy to clean
- High density foam base provides a firm, supportive foundation
- Breathable top cover moves with the patient to reduce shearing
- Soft foam top layer enhances patient comfort
- Zippered cover is removable for easy cleaning



Specifications:

Item	UOM	overall Width		overall Length		overall Thickness		Weight Capacity EVENLY DISTRIBUTED		Bolsters	Soft Heel Zone	Pressure Reduction	Meets 16 CFR 1633	Made in USA	Mattress in months warranty	Cover in months
		inch	cm	inch	cm	inch	cm	lb	Kg							
GEN3980A	1/ea	39	99.06	80	203.20	7	17.78	800	362.87	■	■	■	■	■	60	12
GEN3980B	1/ea	39	99.06	80	203.20	7	17.78	1000	453.59	■	■	■	■	■	60	12
GEN4280A	1/ea	42	106.68	80	203.20	7	17.78	800	362.87	■	■	■	■	■	60	12
GEN4280B	1/ea	42	106.68	80	203.20	7	17.78	1000	453.59	■	■	■	■	■	60	12
GEN4880A	1/ea	48	121.92	80	203.20	7	17.78	800	362.87	■	■	■	■	■	60	12
GEN4880B	1/ea	48	121.92	80	203.20	7	17.78	1000	453.59	■	■	■	■	■	60	12
GEN5480A	1/ea	54	137.16	80	203.20	7	17.78	800	362.87	■	■	■	■	■	60	12
GEN5480B	1/ea	54	137.16	80	203.20	7	17.78	1000	453.59	■	■	■	■	■	60	12
GEN4884H	1/ea	48	121.92	84	213.36	7	17.78	650	294.84	■	■	■	■	■	60	12
650-1501-1633N	1/ea	42	106.68	80	203.20	7	17.78	600	272.16	—	■	■	■	■	60	12
650-1501-4284-1633	1/ea	42	106.68	84	213.36	7	17.78	600	272.16	—	■	■	■	■	60	12

■ Standard
 ■ Optional
 ■ Available
 ■ Not Available
 — Does Not Apply
 🌐 Active Hyperlink



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.