



PATIENT LIFTS AND SLINGS

# MAXICARE BARIATRIC PATIENT LIFTING SYSTEM

When transferring a bariatric patient, safety is key for both the patient and caregiver. The Gendron Patient Lift is designed to offer a safer and robust solution to patient transfers.

**Features**

- Heavy-gauge steel construction
- Four-Point spreader bar with rotation tilts to enhance comfort and safety
- The hand pendant controls all functions of the lift including extending and retracting the legs and the opening and closing of the base
- Lifting and lowering along with the left and right rotational movement are also hand pendant operated
- Lower base height maximizes stability and fits under beds
- 24V DC motor allows quiet smooth operation
- Easy access emergency stop button
- Emergency Power Off button available if needed in event of a power loss
- Powder-coated finish
- Optional digital scale available
- Available in 3 options to support 700 lb (317.5 kg), 800 lb (362.8 kg) and 1,000 (453.5 kg) Weight Capacities

**Warranty:**

- Frame and spreader bar five (5) year limited
- Actuator one (1) year limited



**Specifications:**

	overall Width		overall Length		Minimum Height Spreader bar		Maximum Height Spreader bar		Minimum Width Base		Maximum Width Base		Floor Clearance		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
PL700	85	215.9	61 11/16	156.7	40 13/16	103.7	81 15/16	208.12	53	134.62	85	215.90	4.5	11.43	351	159	700	317.5	1 ea
PL800	85	215.9	61 11/16	156.7	43 9/16	110	78 9/16	199.55	53	134.62	85	215.90	4.75	12.07	354	160.5	800	362.8	1 ea
PL1000	85	215.9	61 11/16	156.7	43 9/16	110	78 9/16	199.55	53	134.62	85	215.90	4.75	12.07	354	160.5	1000	453.5	1 ea