

PRODUCT CATALOG



Maxi Rest Series Beds

Maxi Rest features an adjustable sleeping surface with three width and two length settings. A robust steel two-piece split-frame construction designed for ease of delivery and transport. The 5-inch (12.70 cm) Deluxe Caster Set has 2 Total-Lock and 2 Steer-Lock for ease of maneuverability. Standard 3-Function Positioning Hand Control or optional 5-Function Positioning Hand Control. Synchronized Hi/Lo Motors ensure bed remains level, and battery backup is standard for dependability.



GS4842D

- 800 lb (362.87 kg) capacity EVENLY DISTRIBUTED¹
 - Made in USA — 24 month frame warranty
- [Product Sheet](#)



shown with optional plastic head and foot boards,
standard chrome rail set, and standard 3-function hand control



shown with optional plastic head and foot boards,
plastic rail set, scale system and
5-function hand control

accessories
page 3

Shown with
optional
Full Length
Trapeze



GS4054D

- 1000 lb (453.59 kg) capacity EVENLY DISTRIBUTED¹
 - Made in USA — 24 month frame warranty
- [Product Sheet](#)





Maxi Rest Options and Accessories

available at additional cost

- GL928000 38 inch (96.52 cm) Laminate Wood Grain Head and Foot Board Set
- GL929000 39 inch (99.06 cm) HDPE Plastic Head and Foot Board Set
- GL930000 48 inch (121.92 cm) HDPE Plastic Head and Foot Board Set
- GL927100 Deluxe Plastic Rail Set (LH/RH) [4054 ONLY]
- GL927000 Standard Chrome Rail Set [4842 ONLY]
- GL927020 Standard Chrome Rails w/ 39, 42 & 48 inch Pull-outs
- GL93100 80 inch (203.20 cm) Frame Mount Full Length Trapeze [FITS 80 INCH ONLY]
- GL93180 88 inch (223.52) Frame Mount Full Length Trapeze [FITS 88 INCH ONLY]
- GL93200 5-Function Hand Control (Trend/Rev Trend/Chair)
- GL93300 Full Width Wall Bumper
- GL93600L Egress (Assist) Bar and Left Hand Mounting Pullout
- GL93600R Egress (Assist) Bar and Right Hand Mounting Pullout
- GL94000 8 inch (20.32 cm) Bolster Mattress and Extension Pan Kit
- GS4842DW Bed with Factory Installed Scale System [4842 ONLY]
- GS4054DW Bed with Factory Installed Scale System [4054 ONLY]


GL929000

39" HDPE PLASTIC HEAD AND FOOT BOARD SET



GL927100

DELUXE PLASTIC RAIL SET SHOWN WITH HAND CONTROL



GL927000

STANDARD CHROME RAIL SET SHOWN WITH HAND CONTROL



GS4842DW AND GS4054DW

FACTORY INSTALLED SCALE SYSTEM BEDS




GL93100
FRAME MOUNT FULL LENGTH TRAPEZE
80 inch (203.20 cm) version shown

GL93200

5-FUNCTION TREND/REV TREND HAND CONTROL

Homecare Series Beds

Designed for in-home use and long-term care. Two-piece frame easily disassembles for removing from service and transport. The Gendron Homecare Bed also includes head end side rails. Full-electric with standard 3-Function Positioning hand control.



shown with standard head end side rails



GS3648

- 600 lb (272.16 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 24 month frame warranty
- Product Sheet

TRC 750 Tilt-Recline Chair

Features separate manual hand lever controls which enable the reclining back and seat tilt to function individually. Seat surface is a metal pan covered with 3-inch fabric covered high density foam and is removable if needed. Front riggings are available in Footrests or Full Width Padded Leg Support and Full Width Footboard. Removable Full Length Armrests are standard.



TRC28750 - INCLUDES TILT RECLINER WITH SWING-AWAY FOOTREST, 3 INCH (7.62 CM) PADDED SEAT & BACKREST

TRC28750-FB - INCLUDES TILT RECLINER WITH SWING-AWAY FOOTREST, 3 INCH (7.62 CM) PADDED SEAT & BACKREST, FULL WIDTH PADDED LEG SUPPORT & FOOTBOARD

- 750 lb (340.19 kg) capacity ^{EVENLY DISTRIBUTED}¹
- Made in USA — 36 month warranty

Product Sheet

shown with optional full width foot board (Item No. G27225)

Accessories and Replacement Parts

- G27209 High Density Foam Seat, Black Fabric, 4" High
- G27208 High Density Foam Backrest, Black Fabric 3" Depth
- G27220 Full Width Padded Leg Support
- G27225 Full Width Footboard
- G105605 Headrest Positioning Wings
- G10033-1 Swingaway Footrests-Pair



pneumatic cylinder hand lever control shown



bottom chairs shown with swingaway footrests



Clinical Care Recliner

Features extra-wide seating surface swivel and locking casters for maneuverability. High density foam back with a 42 degree recline and Elevating Legrest support. Recline mechanism located at the rear of the recliner for easy caregiver access.



7150 - LARGE SEAT WIDTH

7155 - MEDIUM SEAT WIDTH

- 850 lb (385.55 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 36 month warranty

 Product Sheet



715x Series Colors

Item	Color
715xP1117	Grey
715xP1199	Dark Blue
715xP180	Brown
715xP183	Black
715xP186	Port Wine
715xP195	Teal

Patient Slings

GTxXL Series

Features washable, durable nylon fabric with four-point attachment. Designed to meet the requirements of Medicare code: E0621



GT4XL

- 800 lb (385.55 kg) capacity EVENLY DISTRIBUTED¹

GT5XL

- 1000 lb (453.59 kg) capacity EVENLY DISTRIBUTED¹

GT5XL - COMMODE OPENING

- 1000 lb (453.59 kg) capacity EVENLY DISTRIBUTED¹

- Made in USA — 6 month warranty

[Web Page](#)

Sling Specifications:

Item	Body Width		Body Length		Weight ² Capacity EVENLY DISTRIBUTED		Size	Accent Color	Commode Opening
	inch	cm	inch	cm	lb	kg			
GT4XL	47.00	119.38	67.00	170.18	800	362.87	M	green	—
GT5XL	47.00	119.38	72.50	184.15	1000	453.59	L	orange	—
GT5XLC	47.00	119.38	72.50	184.15	1000	453.59	XXL	red	■

Maxicare Patient Lift

PLx00 Series

Features hand pendant controls base leg extension and widening for stability as well as patient lifting, transferring left or right, and lowering. Each model is available with optional scales. Designed to meet the requirements of Medicare code: E0635.



PL700

- 700 lb (317.51 kg) capacity EVENLY DISTRIBUTED¹

PL800

- 800 lb (385.55 kg) capacity EVENLY DISTRIBUTED¹

PL1000

- 1000 lb (453.59 kg) capacity EVENLY DISTRIBUTED¹

- Made in USA — 12 month actuator warranty
- 60 month frame/spreader warranty

[Scale Installation Instructions](#)

[Web Page](#)



Lift Specifications:

Item	Base Width ext. legs open minimum		Base Width ext. legs open maximum		Overall Length		spreader bar Height minimum		spreader bar Height maximum		clearance Height base - floor		Weight Capacity EVENLY DISTRIBUTED		Links
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	
PL700	53.00	134.62	85.00	215.90	61.69	156.69	40.81	103.66	81.94	208.13	4.50	11.43	700	317.51	User Manual
PL800	53.00	134.62	85.00	215.90	61.69	156.69	43.31	110.01	78.56	199.54	4.75	12.07	800	362.87	User Manual
PL1000	53.00	134.62	85.00	215.90	61.69	156.69	43.31	110.01	78.56	199.54	4.75	12.07	1000	453.59	User Manual





Portable Patient Helper Trapeze

Helps position patient in bed and transfer from bed to other equipment with minimal assistance. Features durable robust steel frame with casters and floor stand base that folds for transportability.



GH935850

- 850 lb (385.55 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 36 month frame warranty

[Product Sheet](#)

Shower Commode

Features a full width push bar and rust-resistant, powder-coated stainless-steel frame. Available with footrests.



5226 - 26 INCH SEAT WIDTH

5228 - 30 INCH SEAT WIDTH

- 750 lb (340.19 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 60 month frame and 6 month parts warranty

 Product Sheet



Stationary Commode

Features sturdy powder-coated steel frame with adjustable, removable armrest for lateral transfer.



523326 - 26 INCH SEAT WIDTH

523330 - 30 INCH SEAT WIDTH

- 750 lb (340.19 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 60 month frame and 6 month parts warranty

 Product Sheet



Transport Stretchers

Feature two-position vertical side rails with side-mounted lever brakes equipped with steer/brake functions. Four IV sockets with IV pole storage under stretcher top for optional telescoping IV pole Item No. 10122. Includes pad, side rails, and fowler backrest.

Heavy Duty 1190 Series

Has steel powder-coated painted construction with foot pedal controlled hydraulic Fowler backrest. Two-position fold-down side rails. High-density pad with nylon cover provides durability and patient comfort.



G1190

- 1000 lb (453.59 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 24 month warranty

Product Sheet



General Duty 1050 Series

Carbon steel with carbon steel sheet metal construction and articulation for head section provides safe transport. A high-density litter pad with low-shear fabric ensures durability.



G1050

- 350 lb (158.76 kg) capacity EVENLY DISTRIBUTED¹
- Made in USA — 24 month warranty

Product Sheet



Regency Bariatric Wheelchairs

450 and 525 Series

Are designed for optimal bariatric patient transport and built with a fortified, heavy-duty frame and reinforced at all critical load-bearing points. Multiple seat widths (24"-30") enables institutions and caregivers to accommodate a variety of patients. Comfortable, easy-to-clean black nylon upholstery. Removable desk-length or full-length arms with padded armrests. Features black handrims and non-marking polyurethane tires.



Standard Fixed-Back Series

Adjustable-length swingaway footrests with durable black footplates. Designed to meet the requirements of Medicare code: K0007.



450 SERIES

- 450 lb (204.12 kg) capacity EVENLY DISTRIBUTED¹

525 SERIES

- 525 lb (238.14 kg) capacity EVENLY DISTRIBUTED¹

- Made in USA — 60 month warranty

[Product Sheet](#)



*shown: fixed back with removable full arms
- other configurations available*

Reclining Bariatric Series

Dual pneumatic cylinders facilitate full reclining back of this robust, steel reinforced frame. Removable full support headrest and arms. Anti-tippers and elevating leg rests are standard. Designed to meet the requirements of Medicare code: E1226.



*shown with removable desk arms
- other configurations available*



Available Configurations:

Series	Regency	Meets HCPCS	24 inch (60.96 cm)	26 inch (66.04 cm)	28 inch (71.12 cm)	30 inch (76.20 cm)	15.5 inch (39.37 cm)	17.5 inch (44.45 cm)	Removable Full Arms	Removable Desk	Adjust. Standard Foot Rests	Elevating Leg	Removable Head	Anti-Tippers
450xxxx	Fixed Back	K0007	■	■	■	■	—	■	—	■	■	—	—	■
450Rxxxx	Reclining	E1226	■	■	■	■	■	■	—	■	■	■	■	■
525xxxx	Fixed Back	K0007	■	■	■	■	■	■	■	■	■	—	—	■
525Rxxxx	Reclining	E1226	■	■	■	■	■	■	—	■	■	■	■	■

■ Standard

■ Available

■ Optional

— Not Applicable/Available



Regency Bariatric Fixed Back Wheelchair

5600 Series

Designed for optimal bariatric comfort and built with a reinforced, heavy-duty frame, and is reinforced at all critical load-bearing points. Available with adjustable-length swingaway footrests or elevating legrests with durable black footplates.



5600 SERIES

- 600 lb (272.16 kg) capacity ^{EVENLY DISTRIBUTED}¹
- Made in USA — 60 month warranty
- Product Sheet

shown with removable desk arms
- other configurations available



Regency Bariatric Wheelchair

6700 and 6800 Series

Built with dual square cross braces and a heavy-duty frame that is reinforced at all critical load bearing points to support heavier users. Comfortable, easy-to-clean black nylon upholstery. Removable desk- and full-length arm styles with padded armrests.

Standard Fixed-Back Series

Adjustable-length swingaway footrests with durable black footplates. Designed to meet the requirements of Medicare code: K0007.

shown with removable full arms
- other configurations available



Reclining Bariatric Series

Features two pneumatic cylinders which facilitate adjustment from upright to full recline. Removable full support headrest and arms and anti-tippers and elevating leg rests are standard. Designed to meet the requirements of Medicare code: K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K0195 (Elevating Legrests).



shown with removable desk arms
- other configurations available

Available Configurations:

Series	Regency	Meets HCPCS	22 inch (55.88 cm)	24 inch (60.96 cm)	26 inch (66.04 cm)	28 inch (71.12 cm)	30 inch (76.20 cm)	32 inch (81.28 cm)	34 inch (86.36 cm)	18 inch (45.72 cm)	20 inch (50.80 cm)	22 inch (55.88 cm)	16 inch (40.64 cm)	17 inch (43.18 cm)	18 inch (45.72 cm)	19.5 inch (49.53 cm)	Removable Full Arms	Removable Desk Arms	Adjust. Swingaway Foot Rests	Elevating Leg Rests	Removable Head Rest	Anti-Tippers
5600xxx	XL2002 Fixed Back	K0007	—	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
6700xxx	XL2000 Fixed Back	K0007	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
6700Rxxx	XL2000 Recliner	multi ²	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
6800xxx	XL2000 Fixed Back	K0007	—	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
6800Rxxx	XL2000 Recliner	multi ²	—	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■

² see above

■ Standard ■ Available ■ Optional — Not Applicable/Available



6700 SERIES

- 700 lb (317.51 kg) capacity ^{EVENLY DISTRIBUTED}¹

6800 SERIES

- 800 lb (362.87 kg) capacity ^{EVENLY DISTRIBUTED}¹

- Made in USA — 60 month warranty

- Product Sheet

grahamfield

MRI Transport & Transfer Stretcher

Heavy-duty, stainless steel non-magnetic stretcher. With a 2-inch (5.08 cm) high density liter pad with low shear fabric, and articulation to raise and lower the head section. The 2-position swingaway side rails keep patient secure during transport which is facilitated using 4 heavy duty, 6 inch (15.24 cm) locking casters with caster mounted pedal brakes. Stretcher features IV sockets on each corner to use removable, telescopic IV pole.



1000MR

- 350 lb (158.76 kg) capacity ^{EVENLY DISTRIBUTED}¹
- Made in USA — 12 month warranty

[Product Sheet](#)

[Bariatric MRI Platforms Product Sheet](#)



Special Notice: This product (1000MR) is intended for use inside the MRI environment (e.g., in the MR room, close to the scanner) is MR conditional and can be used in the MRI environment according to the following conditions:



- It should not be utilized directly inside the MR system (e.g., inside the bore of the scanner)
- Static magnetic field of 3 - Tesla or less
- Spatial gradient magnetic field of 720 - Gauss/cm or less

MRI Transport Wheelchairs

Constructed of non-magnetic materials and are designed specifically for use in the MRI suite. They should not be utilized directly inside of the MR system. Designed with a distinctive blue color (4650 and 4850 models only) to help staff immediately identify MR safe wheelchair equipment. Features swingaway footrests as well as removable full length armrests for comfort and convenience. Ideal for transferring patients from one department to another.



4000MRQ2

- 350 lb (158.76 kg) capacity ^{EVENLY DISTRIBUTED}¹

4650MR

- 650 lb (294.84 kg) capacity ^{EVENLY DISTRIBUTED}¹

4850MR

- 850 lb (385.55 kg) capacity ^{EVENLY DISTRIBUTED}¹

- Made in USA — 12 month warranty

[Product Sheet](#)

[Bariatric MRI Platforms Product Sheet](#)



Special Notice: Non-clinical testing demonstrated that each product (4650MR / 4850MR) is MR conditional and can be used in the MRI environment according to the following conditions:



- Static magnetic field of 3 - Tesla or less
- Spatial gradient magnetic field of 720 - Gauss/cm or less

¹ refer to user manual

GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

• 1 770.368.4700

Information contained herein is subject to change without notice.

The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field

GF2100047 Rev/A22