



## BARIATRIC PATIENT MOBILITY

# REGENCY 450-525 STANDARD-RECLINING BARIATRIC WHEELCHAIR

The Regency 450-525 Standard Fixed Back and Reclining Bariatric Wheelchairs are manufactured with the same commitment to quality Gendron is known for.

Designed for optimal bariatric patient transport as well as hospital and clinical use, Regency 450 and 525 wheelchairs are built with a reinforced, heavy-duty frame and reinforced at all critical load-bearing points to support up to 450 lb (450 models) or 525 lb (525 models) users. With its multiple seat widths ranging from 24", 26", 28" to 30", with 20" seat depth, and seat-height settings of 15 1/2" (525 models) or 17 1/2" (450 and 525 models), this heavy-duty wheelchair enables institutions and caregivers to accommodate a large variety of patients.

The reclining wheelchair series features a pneumatic-assist reclining back to facilitate upright-to-recline positioning. Rear anti-tip brackets, elevating legrests, and removable headrest are standard.

### Features

- Available in 24", 26", 28" or 30" seat width
- Available in 15 1/2" (525 models) or 17 1/2" (450 and 525 models) seat height
- Standard 20" seat depth
- Robust, steel reinforced frame
- Comfortable, easy-to-clean black nylon upholstery
- Removable desk-length or full-length arms with padded armrests
- Adjustable-length swingaway footrests with durable black footplates for standard models and elevating legrests for reclining models
- Rear 24" x 1" wheels with black handrims and non-marking polyurethane tires
- 450 models: Maximum Weight Capacity 450 lb - EVENLY DISTRIBUTED
- 525 models: Maximum Weight Capacity 525 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code: K0007 for standard Bariatric and K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K0195 (Elevating Legrests)

### Reclining Model Only

- Features pneumatic-assist reclining back to facilitate full recline positioning
- Adjustable-length elevating legrests
- Removable full support headrest and arms

### Warranty

- Five (5) year limited warranty

### Specifications

Wheelchair Model	Patient Weight Capacity EVENLY DISTRIBUTED		overall Width		folded Width		seat width		overall depth		seat depth		overall height		seat height		back height		weight without front rigging	
	lb	kg	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg
450 Standard Fixed Back	450	204.1	33, 35,	83.8, 88.9,	17	43.2	24,	61.0,	43.5	110.5	20	50.8	37	94.0	17.5	44.5	18.5	47.0	63	28.6
450 Recliner			37, 39	94.0, 99.1			26,	66.0,	52.5	133.4			36	91.4			94	42.6		
525 Standard Fixed Back	525	238.1	34, 36,	86.4, 91.4,	19	48.3	28,	71.1,	43.5	110.5	20	50.8	37	94.0	15.5,	39.4,	18.5	47.0	63	28.6
525 Recliner			38, 40	96.5, 101.6			30	76.2	52.5	133.4			17.5	44.5			36	91.4	94	42.6



Regency Standard Bariatric Wheelchair with swingaway footrests and removable full-length arms



Anti-tippers (standard on Recliners)



Regency Reclining Bariatric Wheelchair with elevating legrests and removable desk-length arms



Push handle and brakes (standard on Recliners)

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com) / [www.gendroninc.com](http://www.gendroninc.com)

