



BARIATRIC PATIENT SEATING

BARIATRIC CLINICAL CARE RECLINERS

Designed for patients weighing up to 850 pounds (385.5 kg) when evenly distributed, the Gendron Bariatric Clinical Care Recliner features swivel and locking casters for maneuverability. High density foam Back features a 42 degree recline and Elevating Legrest support. The generous extra-wide seating surface measures 29.25" (74.3cm) wide (Model 7150) or 25.25" (64.1cm) wide (Model 7155) and 19.25" (48.9cm) deep. The recline mechanism is located at the rear of the Recliner for easy caregiver access.



Available Colors

- P1117 - Grey
- P1199 - Dark Blue
- P180 - Brown
- P183 - Black
- P186 - Port Wine
- P195 - Teal

Areas of Use

- Dialysis
- Oncology
- IV Therapy
- Patient Room
- Bariatric Care



Specifications:

	Seat Width		overall Seat Width		Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	lb	Kg	
7150 (color code)	29.25	74.3	38.25	97.2	850	385.5	1/ea
7155 (color code)	25.25	64.1	34.25	87	850	385.5	1/ea

WARNINGS: Cancer and Reproductive Harm - www.P66Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



GF Health Products, Inc. is an
ISO 13485:2016
MDSAP Certified Company
MDSAP
MEDICAL DEVICE SINGLE AUDIT PROGRAM



+1.770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com