



BARIATRIC PATIENT SEATING

BARIATRIC TILT-RECLINE SEATING

The TRC 750 is the choice for tilt and recline positioning of Bariatric patients. Separate manual hand lever controls enable the reclining back and seat tilt to function individually. The seat surface is a metal pan covered with four-inch fabric covered high density foam and is removable if needed. Front riggings are available in Footrests, Elevating Legrests or an optional Full Width Padded Leg Support and Full Width Footboard. Removable Full Length Armrests are standard. Designed for patients weighing up to 750 pounds (340.2 kg).



Back in up position and seat in reverse tilt position.



Reclining back in full upright position and seat in forward tilt position.



Back in full recline position and seat in reverse tilt position.

Accessories

- 100561: High Density Foam Seat, Grey Fabric, 4" High
- 100562: Gel-Foam Pack Insert Seat Cushion, Black, 4" High
- 1002217: Full Width Padded Leg Support
- 100252: Full Width Footboard
- 105605: Headrest Positioning Wings
- 100334: Swing Away Footrests-Pair
- 100361P183: Swing Away Elevating Legrests-Pair

Areas of Use

- Long Term Care
- Respiratory Therapy
- Position Change Therapy
- Patient Room
- Bariatric Care



Push handle and manual control levers. Two levers to control back recline and two levers to control seat tilt.

Back upholstery features adjustable straps to allow upholstery to conform to user

Specifications:

	Seat Width		Seat Depth		Seat pan to floor		Back Height from seat		Back Angle upright/full incline	forward Seat Tilt	reverse Seat Tilt	Weight Capacity EVENLY DISTRIBUTED		UOM
	inch	cm	inch	cm	inch	cm	inch	cm	degrees	degrees	degrees	lb	Kg	
TRC28750	28	71.1	20	50.8	16.25	41.3	35	88.9	10 - 36	5	22	750	340.2	1/ea

WARNING: Cancer and Reproductive Harm - www.P66Warnings.ca.gov

Color variations may occur due to monitor or printer settings.

Note: All images above are shown with optional Full Width Footboard and Full Width Padded Leg Support



GF Health Products, Inc. is an ISO 13485:2016 MDSAP Certified Company



© 2020, GF Health Products, Inc. All Rights Reserved. Gendron is a registered trademark of GF Health Products, Inc.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.grahamfield.com



GF2000247RevA30