



## BARIATRIC PATIENT MOBILITY

### REGENCY 6700-6800 STANDARD-RECLINING BARIATRIC WHEELCHAIR

The Regency 6700-6800 Standard and Reclining Bariatric Wheelchair is manufactured with the same commitment to quality Gendron is known for.

Designed for optimal bariatric patient transport as well as hospital and clinical use, Regency 6700 and 6800 wheelchairs are built with a dual square crossbrace heavy-duty frame and reinforced at all critical load-bearing points to support up to 700 lb (6700 models) or 850 lb (6800 models) users. With its multiple seat widths ranging from 22", 24", 26", 28", 30", 32" to 34", most with 18", 20" or 22" seat depths, and low-seat-height settings of 16", 17" or 18" (factory-set at 18"), this heavy-duty wheelchair enables institutions and caregivers to accommodate a large variety of patients.

The reclining wheelchair series features a pneumatic-assist reclining back to facilitate upright-to-recline positioning. Rear anti-tip brackets and removable headrest and arms are standard.

#### Features

- Available in 22", 24", 26", 28", 30", 32" or 34" seat width
- Available in 16", 17" or 18" seat height
- Available in 18" 20" or 22" (6700 models) and 20" or 22" (6800 models) seat depth
- Robust, steel reinforced frame
- Comfortable, easy-to-clean black nylon upholstery
- Removable desk- and full-length arm styles with padded armrests
- Adjustable standard footrests with durable black footplates and elevating legrests for reclining wheelchair models
- Rear 24" x 1" wheels with black handrims and non-marking polyurethane tires
- 6700 models: Maximum Weight Capacity 700 lb - EVENLY DISTRIBUTED
- 6800 models: Maximum Weight Capacity 850 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code: K0007 for standard Bariatric and K0007/E1226 (Recliner Back) + E0966 (Head Rest Extension) + E0971 (Anti-tipping device) + K0195 (Elevating Legrests)

#### Reclining Model Only

- Features pneumatic-assist reclining back to facilitate full recline positioning
- Removable full support headrest and arms
- Rear anti-tippers standard

#### Warranty

- Five (5) year limited warranty

#### Specifications

Wheelchair Model	Patient Weight Capacity EVENLY DISTRIBUTED		overall Width		folded Width		seat width		overall depth		seat depth		overall height		seat height**		back height		weight without front rigging	
	lb	kg	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg
6700 Fixed Back	700	317.5	31, 33,	78.7, 83.8,	22	55.9	22, 24,	55.9, 61.0,	43	109.2	18,	45.7,	35	88.9	16,	40.6,	18	45.7	72	32.7
6700 Recliner	700	317.5	35, 37,	88.9, 94.0,	22	55.9	26, 28,	66.0, 71.1,	53	134.6	20,	50.8,	53	134.6	17,	43.2,	18	45.7	101	45.8
6800 Fixed Back	850	385.6	39, 41,	99.1, 104.1,	22	55.9	30, 32,	76.2, 81.3,	43	109.2	20,	50.8,	35	88.9	16,	40.6,	18	45.7	85	38.6
6800 Recliner	850	385.6	43	109.2	22	55.9	34	86.4	53	134.6	22	55.9	53	134.6	17,	43.2,	18	45.7	105	47.6

\* If seat depth is 22" (55.9 cm), then overall depth is 46" (116.8 cm)

\*\* Seat height is factory-set at 18" (45.7 cm)



Regency Standard Bariatric Wheelchair with footrests and removable full-length arms



Dual square crossbrace and heavy-duty frame (standard)



Anti-tippers (standard on Recliners)



Regency Reclining Bariatric Wheelchair with elevating legrests and removable desk-length arms



Push handle and brakes (standard on Recliners)

WARNING: Cancer and Reproductive Harm - www.P65Warnings.ca.gov

Color variations may occur due to monitor or printer settings.



+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. [www.grahamfield.com](http://www.grahamfield.com) / [www.gendroninc.com](http://www.gendroninc.com)

