



SUPPORT SURFACES

BARIATRIC SKIN PROTECTION CUSHIONS

The Bariatric Skin Protection Cushions provides pressure distribution and weight equalization to meet unique demands of the bariatric patient.

Available with Foam w/Gel Layers:

- 4" cushion
- Medium density foam base
- Non-aqueous specialized visco-elastic gel core layer reduces tissue interface pressure at critical points
- Proprietary gel is a flexible solid that will not leak
- Gel is placed at three pressure zones
- Visco-elastic memory foam topper for superior comfort and pressure relief
- Positions the user with a moderate slope from front to back

Available with 3 Layer Foam:

- 4" cushion
- High density foam base for structure
- High density foam middle layer for support
- Medium density top foam layer for comfort



Cover Features

- Fluid-resistant stretch top cover, reduces shear and friction
- Antimicrobial
- Flame retardant
- Removable and washable
- Non-skid base
- Strap to secure cushion firmly to the chair

Specifications

Item	UOM	Foam w/ Gel Layer	3 Layer Foam Construction	Meets HCPCS	Overall Width		Overall Length		Overall Thickness		Weight Capacity EVENLY DISTRIBUTED		Made in USA	Warranty in months
					inch	cm	inch	cm	inch	cm	lb	kg		
20252218B	1/ea	■		E2604	22	55.9	18	45.7	4	10.2	650	294.8	■	18
20252418B	1/ea	■		E2604	24	61	18	45.7	4	10.2	650	294.8	■	18
20252619B	1/ea	■		E2604	26	66	19	48.3	4	10.2	650	294.8	■	18
20253019B	1/ea	■		E2604	30	76.2	19	48.3	4	10.2	650	294.8	■	18
G2418-750	1/ea		■		24	61	18	45.7	4	10.2	750	340.19	■	18
G2618-750	1/ea		■		26	66	18	45.7	4	10.2	750	340.19	■	18
G2420-750	1/ea		■		24	61	20	50.8	4	10.2	750	340.19	■	18
G2620-750	1/ea		■		26	66	20	50.8	4	10.2	750	340.19	■	18
G2820-750	1/ea		■		28	71.1	20	50.8	4	10.2	750	340.19	■	18
G3020-750	1/ea		■		30	76.2	20	50.8	4	10.2	750	340.19	■	18



Images are gel foam model



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2008, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field

GF-0800128Rev.024