



MAXI REST BARIATRIC BED

MAXI REST 4054 SERIES

Features:

- Robust steel two-piece split-frame construction for ease of delivery and transport
- Bed deck width is expandable to 39", 48", or 54"
- Bed deck length is adjustable to 80" or 88"
- 5" Deluxe Caster Set includes 2 Total-Lock and 2 Steer-Lock for ease of maneuverability
- Standard 3-Function Positioning hand control
- Optional 5-Function Positioning hand control (Trendelenburg / Reverse Trendelenburg and Chair Positioning)
- Synchronized Hi/Lo Motors ensure bed remains level
- Battery backup standard for dependability
- Any height rolling position for ease of use
- Optional Intuitive side rail system for both head and foot sections
- Optional laminate wood grain or plastic molded head and foot board options
- Optional factory-installed weight scale system
- 1,000 lb (453.59 kg[†]) weight capacity¹ –EVENLY DISTRIBUTED

4045 Maxi Rest Bed Product Info



Options:	Maxi Rest Bed Only	Split Frame	Rail Mounting at head end	Scale
GS4054D	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
GS4054DW	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

★ options available at additional cost

Specifications:

Series	Bed Deck Width Positions						Bed Deck Length Positions				Bed Deck Height Range				Weight ¹ Capacity EVENLY DISTRIBUTED		Standard	Optional	Motors	Battery Backup
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg [†]	functions	features	emergency	
4054	39.00	99.06	48.00	121.92	54.00	137.16	80.00	203.20	88.00	223.52	13.50	34.29	27.50	69.85	1,000	453.59	3	5	4	<input checked="" type="checkbox"/>

[†] standard rounding method used
 ¹ maximum weight capacity - evenly distributed - refer to user manual
 ⁸ approximate measurement



MAXI REST 4054 SERIES



Maxi Rest Options and Accessories

sold separately at additional cost

- GL928000 38 inch (96.52 cm) Laminate Wood Grain Head and Foot Board Set
- GL929000 39 inch (99.06 cm) HDPE Plastic Head and Foot Board Set
- GL930000 48 inch (121.92 cm) HDPE Plastic Head and Foot Board Set
- GL927100 Deluxe Plastic Rail Set (LH/RH)
- GL927000 Standard Painted Head End Rails
- GL927020 Standard Painted Foot End Rails w/ 39, 42 & 48 inch Pull-outs
- GL93100 80 inch (203.20 cm) Frame Mount Full Length Trapeze [FITS 80 INCH ONLY]
- GL93180 88 inch (223.52 cm) Frame Mount Full Length Trapeze [FITS 88 INCH ONLY]
- GL93200 5-Function Hand Control (Trendelenburg/Rev Trendelenburg/Chair Positioning)
- GL93300 Full Width Wall Bumper
- GL93600L Egress (Assist) Bar and Left Hand Mounting Pullout
- GL93600R Egress (Assist) Bar and Right Hand Mounting Pullout
- GL94000 8 inch (20.32 cm) Bolster Mattress and Extension Pan Kit


GL929000
39 INCH HDPE PLASTIC HEAD AND FOOT BOARD SET



GL927000
STANDARD PAINTED RAIL SET SHOWN WITH HAND CONTROL



GL927100
DELUXE PLASTIC RAIL SET SHOWN WITH HAND CONTROL




Warranty

- Main Frame and welds:.....Five years
- Control Box and Actuator Motors:.....Two years
- Hand Control Pendant, Staff Control, and Cabling:.....Two years
- Head Board and Foot Board:.....One year
- All other durable components not listed above:.....Two years

GL93200
5-FUNCTION TREND/REV
TREND HAND CONTROL



GL93100
FRAME MOUNT FULL LENGTH TRAPEZE
80 inch (203.20 cm) version shown



GS4054DW
BED WITH FACTORY INSTALLED SCALE SYSTEM



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2022, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.

graham
field

GF200023Rev/C2.4