




PERSONAL CARE AND AMBULATORY



# SHOWER COMMODE CHAIR

**Features**

- 750 lb maximum weight capacity (Evenly Distributed)
- Rust-resistant, powder-coated stainless steel frame
- Removable armrest
- 2 locking rear casters
- Full-width push bar for easy transport assist
- Tensioning back support straps for better patient positioning
- Mildew- and bacteria-resistant back support fabric
- Pail is removable from the back



**Specifications**

Bariatric Shower Commode Chair Model	Patient Weight Capacity		Overall Width		Seat Width		Seat Depth		Seat-to-Floor Height		Frame Material	Arm Style	Footrest Style
	EVENLY DISTRIBUTED		inch	cm	inch	cm	inch	cm	inch	cm			
	lb	kg											
5226	750	340.2	35.3	89.7	26	66.0	20	50.8	21.7	55.1	Coated Stainless Steel	Removable	Removable
5228	750	340.2	35.3	89.7	28	71.1	20	50.8	21.7	55.1	Coated Stainless Steel	Removable	Removable



GF Health Products, Inc. is an  
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.  
The most current and complete product information can be found on our website.

[www.gendroninc.com](http://www.gendroninc.com) | [www.grahamfield.com](http://www.grahamfield.com)



© 2020, GF Health Products, Inc. All Rights Reserved.  
Gendron is a registered trademark of GF Health Products, Inc.