



SHOWER COMMODE CHAIR

Features

- 750 lb (340.19 kg) weight capacity ^{EVENLY DISTRIBUTED}¹
- Rust-resistant, powder-coated stainless steel frame
- Removable armrest
- 2 locking rear casters
- Full-width push bar for easy transport assist
- Tensioning back support straps for better patient positioning
- Mildew- and bacteria-resistant back support fabric
- Pail is removable from the back

Product Info



Specifications

Item	Overall Width		Overall Depth		Overall Height		Seat Width		Seat Depth		Seat Height		Back Height above seat		Shipping Weight ⁸				Weight ¹ Capacity EVENLY DISTRIBUTED	
	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	without pallet		with pallet		lb	kg†
5226	36.00	91.44	43.00	109.22	40.00	101.60	26.00	66.04	20.00	50.80	20.00	50.80	18.00	45.72	59	26.76	99	44.91	750	340.19
5228							28.00	71.12							60	27.22	100	45.36		

† standard rounding method used ¹ maximum weight capacity - evenly distributed - refer to user manual ⁸ approximate measurement



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2020, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.