



BARIATRIC PATIENT MOBILITY



REGENCY WHEELCHAIR

5600 SERIES FIXED BACK WHEELCHAIRS


Designed for optimal bariatric patient transport as well as hospital and clinical use, the Regency 5600 wheelchairs are built with a heavy-duty frame reinforced at all critical load-bearing points to support a wide variety of patients. With its available seat widths ranging from 24" - 28", this heavy-duty wheelchair accommodates multiple options for patient comfort and mobility.

Features:

- Robust, steel-reinforced frame
- Comfortable, easy-to-clean black nylon upholstery
- Adjustable-length swingaway footrests or elevating legrests with durable black footplates
- Removable desk-length or full-length arms with padded armrests
- Rear 24" x 1" wheels with black handrims and non-marking polyurethane tires
- 600 lb (272.20 kg) Weight Capacity - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code: K0007



5600 fixed-back model with removable desk arms (shown) and adjustable swingaway footrest (shown) — other configurations available



Regency 5600 Bariatric Fixed-Back Wheelchair 5600 Series

Item	Meets HCPCS	Overall Chair Width		Overall Chair Depth		Overall Chair Height		Folded Chair Width		Seat Width		Seat Depth		Seat Height		Back Height		Wheel Diameter		Caster Diameter		Product Weight w/o front rigging		Weight Capacity EVENLY DISTRIBUTED		Detach. Full arms	Detach. Desk	Adjust. Swingaway Foot rests	Elevating Leg	Rear Anti-Tippers
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	lb	kg	lb	kg							
56241820	K0007	33.00	83.82	43.50	110.49	37.00	93.98	17.00	43.18	24.00	60.96	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	■	—	■
56241830	K0007	33.00	83.82	43.50	110.49	37.00	93.98	17.00	43.18	24.00	60.96	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	—	■	■
56241840	K0007	33.00	83.82	43.50	110.49	37.00	93.98	17.00	43.18	24.00	60.96	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	■	—	■
56241850	K0007	33.00	83.82	43.50	110.49	37.00	93.98	17.00	43.18	24.00	60.96	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	—	■	■
56261820	K0007	35.00	88.90	43.50	110.49	37.00	93.98	17.00	43.18	26.00	66.04	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	■	—	■
56261830	K0007	35.00	88.90	43.50	110.49	37.00	93.98	17.00	43.18	26.00	66.04	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	—	■	■
56261840	K0007	35.00	88.90	43.50	110.49	37.00	93.98	17.00	43.18	26.00	66.04	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	■	—	■
56261850	K0007	35.00	88.90	43.50	110.49	37.00	93.98	17.00	43.18	26.00	66.04	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	—	■	■
56281820	K0007	37.00	93.98	43.50	110.49	37.00	93.98	17.00	43.18	28.00	71.12	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	■	—	■
56281830	K0007	37.00	93.98	43.50	110.49	37.00	93.98	17.00	43.18	28.00	71.12	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	—	■	—	■	■
56281840	K0007	37.00	93.98	43.50	110.49	37.00	93.98	17.00	43.18	28.00	71.12	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	■	—	■
56281850	K0007	37.00	93.98	43.50	110.49	37.00	93.98	17.00	43.18	28.00	71.12	18.00	45.72	19.50	49.53	18.50	46.99	24.00	60.96	8.00	20.32	63	28.58	600	272.16	■	—	—	■	■

■ Standard
 ■ Optional
 ■ Available
 ■ Not Available
 ■ Accessory
 — Does Not Apply
 Active Hyperlink

¹ approximate product weight — as seat width increases, so does product weight, however slightly
³ adjustable height full arm accessory may be ordered separately at additional cost - Item No. G10063

[5600 Web Page](#)



GF Health Products, Inc. is an
ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.gendroninc.com | www.grahamfield.com



© 2021, GF Health Products, Inc. All Rights Reserved.
Gendron is a registered trademark of GF Health Products, Inc.