


TRC 750 TILT-RECLINE BARIATRIC CHAIR



Push handle and manual control levers.
Two levers to control back recline and two
levers to control seat tilt.

Specifications

Seat Width: 28"

Seat Depth: 20"

Seat Pan to Floor: 16.25"

Back Height: 35"

Back Angle Upright/Full Incline: 10-36 degrees

Seat Tilt Forward: 5 degrees

Seat Tilt Reverse: 22 degrees

Seat Cushion: Four inch high density foam

Caster: Two rear locking casters

Weight Capacity: 750 pounds

Note: all images shown with optional full width
foot board and full width leg guard



Back upholstery features adjustable straps
to allow upholstery to conform to the user..

GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
Phone: 800.537.2521 www.gendroninc.com