

TRC 750 TILT-RECLINE BARIATRIC CHAIR



New from Gendron, the TRC 750 is an alternate choice for tilt and recline positioning. Separate manual hand lever controls enable the reclining back and seat tilt to function individually. The seat surface is metal pan covered with four-inch fabric covered high density foam and is removable if needed. Front riggings are available in footrests and elevating legrests. An optional full width padded leg support and full width footboard are available. Removable full length arm rests are standard. Dimensions and specifications are shown on the back of this page.



Reclining back in the upright position and seat in reverse tilt position.



Reclining back in the full recline position and seat in reverse tilt position.



Reclining back in full upright position and seat in forward tilt position.

GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
Phone: 800.537.2521 www.gendroninc.com