

Safety Sheets



When bed side rails cannot be used, Safety Sheets are the alternative. The patient/resident will know they are close to the edge of the bed when touching the head and/or foot end bolster.

Gendron Safety Sheets are fully sewn, moisture and air permeable and include tie down straps for secure fit. Foam bolster is concealed in a zipper closure. Designed to fit a mattress or support surface up to 10" thick.

Sizes Available

42"W x 80"L

48"W x 80"L

42"W x 88"L

48"W x 88"L


GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
Phone: 419.636.0848 www.gendroninc.com
Toll free: 800.537.2521 Fax: 419.636.1654

Patient Lift Slings

Patient lift slings designed for the bariatric patient. Four point attachment. All Gendron Slings are washable. Patient weight capacity from 400 to 1,000 pounds.

Dimensions Approx. (inches)	1	2	3
GT4XL (400-800 lbs.)	64.5	52	33.5
GT5XL (400-1,000 lbs.)	71.5	57	34.5
GT5XLC (400-1,000 lbs.)	70	58	35



Wheelchair Cushions

Wheelchair cushions designed for the bariatric patient. Available with or without rigid base. Four inch high density foam for necessary support. Attaching strap included. Moisture resistant zipper nylon cover.

Sizes Available

24"W x 18"L
 24"W x 20"L
 26"W x 18"L
 26"W x 20"L
 28"W x 20"L
 30"W x 20"L

Special sizes available upon request.



GENDRON 
 Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
 Phone: 419.636.0848 www.gendroninc.com
 Toll free: 800.537.2521 Fax: 419.636.1654