


GENDRON®



**TOTAL SOLUTIONS FOR BARIATRIC
PATIENT CARE**


Acute Care Beds



Elite Power Assist Drive Bed

Power Assist Drive acute care bariatric bed reduces the risk of caregiver injury by eliminating manually pushing the bed. Bed deck with adjustment from 39"-48"-54" wide, 1,000 pound patient weight capacity. Five function articulation, CPR for the back rest, head and foot end side rails, four wheel central locking are standard equipment. The bed deck is length adjustable for taller patients. Integral scale system controls located at the foot end. Power Assist Drive controls located at the head end. Shown with optional mattress and patient helper.

Model 5154PDX.



Sentinel II Bariatric Low Bed

Features patient weight capacity of 850 pounds and 9.75" low bed deck to floor height for patients at risk of falling. Five function bed deck articulation, bed deck length adjustment, width adjustable to 48", CPR back rest release, bed exit alarm, hand pendant control and key pad control located at the head end side rails, scale system and care giver controls located at foot end.

Model 3875



Elite Bariatric Bed

Features 1,000 pound patient weight capacity and bed deck width adjustment from 39"-48"-54". The bed deck is length adjustable from 84" to 92" for taller patients. CPR release, scale system, four wheel central locking, five function bed deck articulation, all standard equipment.

Model 5154DX

Maxi Rest 42 Beds

Featuring 39"-42"-48" bed deck width adjustment and 800 pound patient weight capacity to accommodate the bariatric user. Two piece split frame, choice of three or five functions, fully electric with hand pendant control, optional scale system, two head/foot board color choices. Available with or without integral battery back up. The bed is length adjustable from 80" to 88" for taller patients.

Model MC 4842



Bariatric Home Care / Long Term Care Bed

Designed for in-home use and long term care. Two piece frame easily disassembles for removing from service and transport. Sleep surface measures 48" wide and 84" long. Three function pendant control and head end side rails are included. Patient weight capacity is 650 pounds.

Model 3648



Bariatric Low Bed

For the bariatric user at risk for falling. Three function, 850 pound weight capacity, two piece frame for portability, 9.75" bed deck low position, 84" long, 36"-42"-48" bed deck width adjustment and optional head and foot end side rails.

Model 3150





Complete Care Bed

For in-home or long term care, features 42" wide by 84" long sleep surface, 650 pound weight capacity, two side rail choices and pendant control. Optional foot end controls (as shown) available.

Model 3642



Maxi Rest Bed

To meet most all needs for the bariatric user. Two piece or one piece frame available. Patient weight capacity 750 pounds or 1000 pounds. Three or five function articulation, with or without battery back up.

Head end side rails, head and foot boards are standard. Bed Deck width adjustment 39"-48" (for 750 pound rated beds) or 39"-48"-54" for beds rated at 1,000 pounds. Hand pendant control.

Model MC4748 Weight capacity 750 pounds

Model MC4054 Weight capacity 1,000 pounds

All beds shown are designed and manufactured by Gendron. All Gendron beds are designed to meet entrapment requirements for patient safety. We're here to help you choose the correct bed for the individual. Please contact us with any questions or specific special requirement.

For detailed information and diagrams showing entrapment zones, visit www.gendroninc.com



Back in up position and seat in reverse tilt position.



Back in up position and seat in forward tilt position.

The TRC 750 Chair

An alternative choice for tilt and recline positioning. Separate manual hand levers individually control pneumatic cylinders that enable the reclining back and seat tilt to function independently. The metal seat pan surface will accommodate a seat cushion measuring 28" wide by 20" deep. Front riggings are available in footrests or elevating leg rests. Optional full width foot board and padded leg support, 4" foam seat cushion and 4" foam/gel cushion are available.

Specifications:

28" seat width

20" seat depth

Seat pan to floor: 16.25"

Back angle: upright / full recline 10-36 degrees

Seat tilt forward: 5 degrees

Seat tilt reverse: 22 degrees

Weight capacity: 750 pounds

Casters: two rear locking



Two sets of levers control back recline and seat tilt, while adjustable straps allow upholstery to conform to the user.



Back in full recline position and seat in reverse tilt position.

Bariatric Patient Room Seating

Portable wheeled bariatric patient recliner chair features reclining back and leg rest that elevates when the back is reclined. The back can be positioned from full upright to fully reclined. Two models available: 29.25" seat width and 25.25" seat width. Both feature 850 pound weight capacity. Recline mechanism is located at the rear of the chair for attendant access.

Contact us for upholstery color choices.

Model 7150 29.25" seat width

Model 7155 25.25" seat width





Maxi Care Bariatric Patient Lift

When transferring the bariatric patient safety is first for the patient and caregiver. The hand pendant controls base leg extension and widening for stability. Lifting, transferring left and right and lowering the patient are also hand pendant controlled. The actuator controlling raising/lowering features a manual release to lower the patient in the event of power loss. Three models feature 700, 800 and 1,000 pound weight capacities. Each is available with optional weighing scale.

Model PL700 700 Pound capacity

Model PL800 800 Pound capacity

Model PL1000 1,000 Pound capacity



Bariatric Walker

Folding bariatric walker feature 750 pound weight capacity.

Front wheels are optional.

Model 5176



Bariatric Portable Patient Helper

When the bed isn't equipped with a trapeze frame to assist with patient movement, choose the patient helper.

Two models, one rated at 850 pounds, one rated at 600 pounds.

Specifications:

35" boom extension

Overall height 72"

Base measures 26" wide, 61" long

Two wheels for portability

Model 5160/850 pound capacity

Model 5165/600 pound capacity



Bariatric Shower Commode Chair

Gendron's new shower commode chair features tensioning back support straps to better position the individual. The back upholstery allows water to pass through and is mildew and bacteria resistant. Arm rests are removable for lateral transfer.

The commode pail is removable from the rear of the chair. Two locking rear casters. The frame is constructed of rust resistant powder coated stainless steel. Full width push bar at the top of the back.

Model	Seat Width	Seat Depth	Seat to Floor Ht.	Arm Style	Overall Width	Footrests	Weight Capacity
5226	26"	20"	21.7"	removable	35.3"	removable	750 Pounds
5228	28"	20"	21.7"	removable	35.3"	removable	750 Pounds

Stationary Bariatric Commode

Durable steel powder coated frame, resists harsh cleaning solutions and corrosion, adjustable arm rests and multiple seat width and weight capacity choices. Pail removes from the rear of the commode. Easily cleaned and sanitized.

Model choices

26" wide seat surface, height adjustable, 750 pound capacity (Model 5233-26)

30" wide seat surface, height adjustable, 750 pound weight capacity (Model 5233-30)

30" seat width, non-height adjustable, 1,000 pound weight capacity, (Model 5233-30 NA)





Regency 450 Bariatric Wheelchairs

Mobility for the individual weighing up to 450 pounds. Four seat width choices and low seat to floor height

24", 26", 28" and 30" seat widths

Seat depth 20"

Seat to floor height 17.5"

Black nylon upholstery



Regency 450 Recliner Wheelchairs

Bariatric recliner wheelchairs feature seat options specifically for the bariatric user. Dual hand release levers to lower and raise the back. Infinite position recliner from up-right to full recline utilizing two pneumatic cylinders to assist with back positioning. Rear anti-tip brackets, removable arms, black nylon upholstery are standard.

24", 26", 28" and 30" seat widths

Seat depth 20"

Seat to floor height 17.5"

Weight capacity 450 pounds

Removable head rest

Black nylon upholstery

Dual pneumatic controls

The Regency 525 family of bariatric wheelchairs.

The Regency 525 family is designed for those individuals needing a deeper seat, more seat width options and lower seat to floor height for proper positioning. The mid-range weight capacity accommodates patients requiring extra seat width and depth. The Regency 525 is available in permanent back style and high back recliner.



Regency 525 Wheelchairs

Permanent back, 18" high

Seat depth of 20" on all models

Available seat widths: 24", 26", 28", 30"

Weight capacity: 525 pounds

Two choices for seat to floor height : 15.5" or 17.5"



Specially designed tensioning back upholstery on each recliner allows back upholstery adjustment to provide the proper fit for the patient. Adjustable straps allow the back upholstery to conform to the individual.



Regency 525 Recliners

Seat depth of 20" on all models

Available seat widths: 24", 26", 28", 30"

Weight capacity: 525 pounds

Two choices for seat to floor height : 15.5" or 17.5"

Regency XL 2000/2002 Wheelchairs

Manual wheelchairs designed for the bariatric patient featuring seat width and depth options, weight capacity choices, permanent back style or recliner wheelchairs.

Regency XL2002

Weight capacity 600 pounds
Seat widths 24" 26" 28"
Seat depth 18"
Seat to floor height: 19.5"
Desk length or full length arm rests
Black nylon upholstery



Regency XL2000 700/750

Weight capacity 700 to 750 pounds
Seat widths 20" 22" 24" 26" 28" 30" 32" or 34"
Seat depth 18" 20" or 22"
Seat to floor height: 16" 17" or 18" (15" w/special rear wheels)
Back Style: permanent, height adjustable 16" 17" or 18"
Cross brace: dual square tube for weight support
Desk length or full length arm rests
Black nylon upholstery



Regency XL2000 850

Weight capacity 850 pounds
Seat widths 26" 28" 30" 32" or 34"
Seat depth 20" or 22"
Seat to floor height: 16" 17" or 18" (15" w/special rear wheels)
Back Style: permanent, 16" 17" or 18"
Desk length or full length arm rests
Black nylon upholstery

Regency XL2000 Recliner Wheelchairs

Weight capacity 700, 750 or 850 pounds
Pneumatic assist back recline to 45 degrees
Cross brace: dual square tube
Back height 36" from seat
Rear anti-tip brackets included
Seat widths: 20" to 34" depending on weight capacity
Seat depths: 18" to 22" depending on weight capacity
Seat to floor height: 16" 17" or 18"
Desk length or full length arm rests
Black nylon upholstery



Patient Transport/Specialty Wheelchairs



Series 1050 General Duty Transport Stretcher

General duty transport stretcher, MR safe
Swing down side rails
Four individual locking casters
Includes pad, side rails and fowler backrest
350 pound weight capacity



Series 1190 Bariatric Transport Stretcher

General duty transport stretcher, MR safe
Swing down side rails
Four individual locking casters
Includes pad, side rails and fowler backrest
1,000 pound weight capacity

MR Safe Patient Transport Equipment



MR Safe Bariatric Wheelchairs

Wheelchairs:

Safe for use in and around the MR suite
Constructed of non-magnetic materials
Full length removable arms, swing away footrests and rear anti tip brackets included
Blue color (for 26" seat width) for easy identification

Model Choices:

20" seat width, 18" deep seat, 350 pound weight capacity, brushed stainless steel frame Model MR4000 Q2
26" seat width, 20" deep, blue frame, 650 pound weight capacity Model 4650 MR
26" seat width, 20" deep, blue frame, 850 pound weight capacity Model 4850 MR

MR Safe test data available upon request



MR Safe Transport Stretcher

General duty transport stretcher, MR safe
Swing down side rails
Four individual locking casters
Includes pad, IV holder, side rails and fowler backrest
350 pound weight capacity



Custom Design Manual and Power Bariatric Wheelchairs

The Regency XLC can be configured to meet the patient's specific needs.

Features include:

Rigid frame power drive or manual design

Power wheelchairs include fully programmable controls and integral electronic wheel locks.

Customize the XLC:

Seat widths from 20" to 40" (please call for widths greater than 40")

Seat depths from 18" to 26"

Seat to floor heights from 14" to 21"

Back height from 16" to 28"

Manual back decline

Power tilt in space

Power elevating leg rests


Seating choices include solid seat pan for adaptive seating, foam padded seating or bariatric pressure relief seat platform

Weight capacity up to 700 pounds (please call for weight capacity greater than 700 pounds)

For configurations and dimensions not listed here, please contact customer service.

Measuring Proper Fit & Wheelchair Construction

(For a detailed measuring guide and product order form visit www.gendroninc.com)



1. Seat distance from floor (Seat Height) _____
 2. Depth of seat from front to rear (Seat Depth) _____
 3. Overall width of seat (Seat Width) _____
 4. Seat to top of back (Back Rest Height) _____
 5. Seat to top of armrest (Armrest Height) _____
- Back Angle (degrees)
If the seat plane is not parallel to the floor. Indicate the distance from the front and rear of seat to floor.
 - Front Edge to Floor
 - Rear Edge to Floor

Please check for environmental access restrictions.
Tell us the maximum overall width, wheel to wheel, to ensure passage through doorways, ramps, etc.

Gendron recommends that specifications to build each wheelchair be supplied by a physician, physical therapist or other qualified medical professional.

CONTACT INFORMATION

Gendron Inc.

520 W Mulberry St.

Bryan OH 43506

419-636-0848

419-636-9261 Fax

800-537-2521 Customer Service

www.gendroninc.com

E Mail: sales@gendroninc.com

