

MAXI REST 42/SR



(Shown here with optional mattress)

MAXI REST 42/SR Features

- NEW side rail system when head and foot end side rails are needed
- Two piece split frame construction for ease of delivery and transport
- Bed deck width adjustment from 39" to 42" to 48"
- Patient weight capacity 800 pounds
- Wood grain or molded head and footboard options
- Bed deck length adjustment from 80" to 88"
- Battery backup included
- Optional scale system



See reverse side for ordering guide.

GENDRON 
Bariatric Rehab Leadership

520 West Mulberry Street Bryan, OH 43506
Phone: 800.537.2521 www.gendroninc.com

MAXI REST 42/SR

MAXI REST 42/SR Ordering Guide

MC 4842 SD-PSR	3 function, wood grain head and foot boards, molded side rails at head end
MC 4842 SD-PSR/2	3 function, wood grain head and foot boards, molded side rails at head and foot end
MC 4842 SD-1-PSR	3 function, molded head and foot boards, molded side rails at the head end
MC 4842 SD-1-PSR/2	3 function, molded head and foot boards, molded side rails at the head and foot end
S (add to bed part number)	Factory installed scale system



520 West Mulberry Street Bryan, OH 43506
Phone: 800.537.2521 www.gendroninc.com